



PIG-X (G-15): sc-100096

BACKGROUND

Several cell surface proteins are attached to the membrane through their C-terminal domain and a glycosylphosphatidylinositol (GPI) moiety. Phosphatidylinositol-glycans (PIGs) are multi-pass transmembrane proteins that localize to the endoplasmic reticulum. PIGs are crucial for the synthesis of very early intermediates in GPI-anchor biosynthesis. PIG-X (Phosphatidylinositol-glycan biosynthesis class X protein) is a 258 amino acid endoplasmic reticular protein that, along with PIG-M, is an essential component of GPI-mannosyltransferase 1, an enzyme that transfers the first of the four mannoses in the GPI-anchor precursors. Due to evidence showing that expression of PIG-M was very low in the absence of coexpressed PIG-X, it is likely that PIG-X stabilizes PIG-M. There are two isoforms of PIG-X that are produced as a result of alternative splicing events.

REFERENCES

1. Yeh, E.T., Kamitani, T. and Chang, H.M. 1994. Biosynthesis and processing of the glycosylphosphatidylinositol anchor in mammalian cells. *Semin. Immunol.* 6: 73-80.
2. Kinoshita, T., Ohishi, K. and Takeda, J. 1997. GPI-anchor synthesis in mammalian cells: genes, their products, and a deficiency. *J. Biochem.* 122: 251-257.
3. Maeda, Y., Watanabe, R., Harris, C.L., Hong, Y., Ohishi, K., Kinoshita, K. and Kinoshita, T. 2001. PIG-M transfers the first mannose to glycosylphosphatidylinositol on the luminal side of the ER. *EMBO J.* 20: 250-261.
4. Kang, J.Y., Hong, Y., Ashida, H., Shishioh, N., Murakami, Y., Morita, Y.S., Maeda, Y. and Kinoshita, T. 2005. PIG-V involved in transferring the second mannose in glycosylphosphatidylinositol. *J. Biol. Chem.* 280: 9489-9497.
5. Ashida, H., Hong, Y., Murakami, Y., Shishioh, N., Sugimoto, N., Kim, Y.U., Maeda, Y. and Kinoshita, T. 2005. Mammalian PIG-X and yeast Pbn1p are the essential components of glycosylphosphatidylinositol-mannosyltransferase I. *Mol. Biol. Cell.* 16: 1439-1448.
6. Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 610276. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. Kim, Y.U., Ashida, H., Mori, K., Maeda, Y., Hong, Y. and Kinoshita, T. 2007. Both mammalian PIG-M and PIG-X are required for growth of GPI14-disrupted yeast. *J. Biochem.* 142: 123-129.
8. Kim, Y.U. and Hong, Y. 2007. Functional analysis of the first mannosyltransferase (PIG-M) involved in glycosylphosphatidylinositol synthesis in *Plasmodium falciparum*. *Mol. Cells.* 24: 294-300.

CHROMOSOMAL LOCATION

Genetic locus: PIGX (human) mapping to 3q29; Pigx (mouse) mapping to 16 B2.

SOURCE

PIG-X (G-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of PIG-X of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-100096 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

PIG-X (G-15) is recommended for detection of PIG-X of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other PIG family members.

Suitable for use as control antibody for PIG-X siRNA (h): sc-78258, PIG-X siRNA (m): sc-152258, PIG-X shRNA Plasmid (h): sc-78258-SH, PIG-X shRNA Plasmid (m): sc-152258-SH, PIG-X shRNA (h) Lentiviral Particles: sc-78258-V and PIG-X shRNA (m) Lentiviral Particles: sc-152258-V.

Molecular Weight of PIG-X: 29 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.