



Dynactin 2 (73-A2): sc-100931

BACKGROUND

Dynactin is a multi-subunit complex that functions as a binding partner for the Dynein microtubule motor. Dynactin-Dynein binding may be required for most, if not all, cytoplasmic Dynein-driven activities and is thought to contribute to the functional diversity of Dynein. Dynactin 2, also known as DCTN2, Dynamitin or DCTN50, is a peripheral membrane protein that is one of many subunits in the Dynactin complex. Like other Dynactin subunits, Dynactin 2 mediates Dynein-organelle binding and helps to regulate chromosome alignment during prometaphase and spindle organization during mitosis. Overexpression of Dynactin 2 disrupts the Dynactin-Dynein complex, thus inhibiting retrograde axonal transport and causing motor neuron degeneration. Additionally, overexpression of Dynactin 2 may disrupt the cell cycle and lead to osteosarcoma, suggesting a possible role for Dynactin 2 in carcinogenesis.

REFERENCES

- Echeverri, C.J., Paschal, B.M., Vaughan, K.T. and Vallee, R.B. 1996. Molecular characterization of the 50 kDa subunit of Dynactin reveals function for the complex in chromosome alignment and spindle organization during mitosis. *J. Cell Biol.* 132: 617-633.
- Berrueta, L., Tirnauer, J.S., Schuyler, S.C., Pellman, D. and Bierer, B.E. 1999. The APC-associated protein EB1 associates with components of the Dynactin complex and cytoplasmic Dynein intermediate chain. *Curr. Biol.* 9: 425-428.
- Merdes, A., Heald, R., Samejima, K., Earnshaw, W.C. and Cleveland, D.W. 2000. Formation of spindle poles by Dynein/Dynactin-dependent transport of NuMA. *J. Cell Biol.* 149: 851-862.
- Karki, S., Tokito, M.K. and Holzbaur, E.L. 2000. A Dynactin subunit with a highly conserved Cysteine-rich motif interacts directly with Arp1. *J. Biol. Chem.* 275: 4834-4839.
- Hoogenraad, C.C., Akhmanova, A., Howell, S.A., Dortland, B.R., De Zeeuw, C.I., Willemsen, R., Visser, P., Grosveld, F. and Galjart, N. 2001. Mammalian Golgi-associated Bicaudal-D2 functions in the Dynein-Dynactin pathway by interacting with these complexes. *EMBO J.* 20: 4041-4054.
- LaMonte, B.H., Wallace, K.E., Holloway, B.A., Shelly, S.S., Ascaño, J., Tokito, M., Van Winkle, T., Howland, D.S. and Holzbaur, E.L. 2002. Disruption of Dynein/Dynactin inhibits axonal transport in motor neurons causing late-onset progressive degeneration. *Neuron* 34: 715-727.
- Uetake, Y., Terada, Y., Matuliene, J. and Kuriyama, R. 2004. Interaction of CEP135 with a p50 Dynactin subunit in mammalian centrosomes. *Cell Motil. Cytoskeleton* 58: 53-66.
- Brocard, C.B., Boucher, K.K., Jedesko, C., Kim, P.K. and Walton, P.A. 2005. Requirement for microtubules and Dynein motors in the earliest stages of peroxisome biogenesis. *Traffic* 6: 386-395.
- Bransfield, K.L., Askham, J.M., Leek, J.P., Robinson, P.A. and Mighell, A.J. 2006. Phenotypic changes associated with Dynactin 2 (DCTN2) over expression characterise SJS-1 osteosarcoma cells. *Mol. Carcinog.* 45: 157-163.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

CHROMOSOMAL LOCATION

Genetic locus: DCTN2 (human) mapping to 12q13.3.

SOURCE

Dynactin 2 (73-A2) is a mouse monoclonal antibody raised against recombinant Dynactin 2 of human origin.

PRODUCT

Each vial contains 100 µg IgG₁ in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Dynactin 2 (73-A2) is recommended for detection of Dynactin 2 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Dynactin 2 siRNA (h): sc-95708, Dynactin 2 shRNA Plasmid (h): sc-95708-SH and Dynactin 2 shRNA (h) Lentiviral Particles: sc-95708-V.

Molecular Weight of Dynactin 2: 50 kDa.

Positive Controls: MCF7 whole cell lysate: sc-2206.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-mouse IgG-FITC: sc-2010 (dilution range: 1:100-1:400) or goat anti-mouse IgG-TR: sc-2781 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.