

AP-2 δ (J-18): sc-101886

BACKGROUND

AP-2 transcription factor family members include AP-2 α , AP-2 β , AP-2 γ , AP-2 δ and AP-2 ϵ , which specifically bind to DNA and regulate transcription of selected genes. AP-2 proteins contain a helix-span-helix motif at their C-terminus and a basic central region that, together, mediate DNA binding and dimerization. AP-2 family members have various roles in apoptosis, development, morphogenesis and cell-cycle control. AP-2 δ , also known as TFAP2D or TFAP2BL1 (transcription factor AP-2 β -like 1), is a nuclear protein and is predominantly expressed in skeletal muscle, brain, small intestine, prostate, placenta and thymus. AP-2 δ binds to DNA as a dimer, associated either as a homodimer or as a heterodimer with other members of the AP-2 family. Distinct from other members of the family, AP-2 δ exhibits a different DNA sequence affinity and lacks the PY motif as well as other critical residues in its transactivation domain. This suggests that AP-2 δ may interact with a separate group of coactivators and transactivate genes differently than the other AP-2 proteins.

REFERENCES

- Zhao, F., Satoda, M., Licht, J.D., Hayashizaki, Y. and Gelb, B.D. 2001. Cloning and characterization of a novel mouse AP-2 transcription factor, AP-2 δ , with unique DNA binding and transactivation properties. *J. Biol. Chem.* 276: 40755-40760.
- Cheng, C., Ying, K., Xu, M., Zhao, W., Zhou, Z., Huang, Y., Wang, W., Xu, J., Zeng, L., Xie, Y. and Mao, Y. 2002. Cloning and characterization of a novel human transcription factor AP-2 β like gene (TFAP2BL1). *Int. J. Biochem. Cell Biol.* 34: 78-86.
- Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 610161. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Zhao, F., Lufkin, T. and Gelb, B.D. 2003. Expression of Tfp2d, the gene encoding the transcription factor Ap-2 δ , during mouse embryogenesis. *Gene Expr. Patterns* 3: 213-217.
- Eckert, D., Buhl, S., Weber, S., Jäger, R. and Schorle, H. 2005. The AP-2 family of transcription factors. *Genome Biol.* 6: 246.
- Wenke, A.K., Rothhammer, T., Moser, M. and Bosserhoff, A.K. 2006. Regulation of integrin α 10 expression in chondrocytes by the transcription factors AP-2 ϵ and Ets-1. *Biochem. Biophys. Res. Commun.* 345: 495-501.
- Kim, J.M., Lee, K.H., Jeon, Y.J., Oh, J.H., Jeong, S.Y., Song, I.S., Kim, J.M., Lee, D.S. and Kim, N.S. 2006. Identification of genes related to Parkinson's disease using expressed sequence tags. *DNA Res.* 13: 275-286.
- Sun, L., Huang, S., Wu, Q., Gu, S., Fu, X., Yu, K., Lu, F., Ji, C., Feng, C., Sun, R., Xie, Y. and Mao, Y. 2007. Identification of genes differentially regulated by transcription factor, AP-2 δ . *Front. Biosci.* 12: 1699-1706.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: TFAP2D (human) mapping to 6p12.3; Tcfap2d (mouse) mapping to 1 A3.

SOURCE

AP-2 δ (J-18) is a purified rabbit polyclonal antibody raised against AP-2 δ of human origin.

PRODUCT

Each vial contains 100 μ g IgG in 1.0 ml PBS with < 0.1% sodium azide, 0.1% gelatin and < 0.02% sucrose.

APPLICATIONS

AP-2 δ (J-18) is recommended for detection of AP-2 δ of mouse, rat, human, dog and zebrafish origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Suitable for use as control antibody for AP-2 δ siRNA (h): sc-95077, AP-2 δ siRNA (m): sc-141133, AP-2 δ shRNA Plasmid (h): sc-95077-SH, AP-2 δ shRNA Plasmid (m): sc-141133-SH, AP-2 δ shRNA (h) Lentiviral Particles: sc-95077-V and AP-2 δ shRNA (m) Lentiviral Particles: sc-141133-V.

Molecular Weight of AP-2 δ : 50 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



AP-2 δ (J-18): sc-101886. Western blot analysis of AP-2 δ expression in fetal muscle tissue extract.

RESEARCH USE

For research use only, not for use in diagnostic procedures.