

AP-2 ϵ (T-25): sc-101887

BACKGROUND

AP-2 transcription factor family members include AP-2 α , AP-2 β , AP-2 γ , AP-2 δ and AP-2 ϵ , which specifically bind to DNA and regulate transcription of selected genes. AP-2 proteins contain a helix-span-helix motif at their C-terminus and a basic central region that, together, mediate DNA binding and dimerization. AP-2 family members have various roles in apoptosis, development, morphogenesis and cell-cycle control. AP-2 ϵ , also known as TFAP2E or AP2E, is a nuclear protein and is predominantly expressed in skin, HeLa cells, primary keratinocytes and immortalized keratinocytes. AP-2 ϵ binds to DNA as a dimer, associated either as a homodimer or as a heterodimer with other members of the AP-2 family. Due to its high level of expression in skin, AP-2 ϵ is believed to play an important role in skin biology.

REFERENCES

- Zhao, F., Satoda, M., Licht, J.D., Hayashizaki, Y. and Gelb, B.D. 2001. Cloning and characterization of a novel mouse AP-2 transcription factor, AP-2 δ , with unique DNA binding and transactivation properties. *J. Biol. Chem.* 276: 40755-40760.
- Tummala, R., Romano, R.A., Fuchs, E. and Sinha, S. 2003. Molecular cloning and characterization of AP-2 ϵ , a fifth member of the AP-2 family. *Gene* 321: 93-102.
- Feng, W. and Williams, T. 2003. Cloning and characterization of the mouse AP-2 ϵ gene: a novel family member expressed in the developing olfactory bulb. *Mol. Cell. Neurosci.* 24: 460-475.
- Wang, H.V., Vaupel, K., Buettner, R., Bosserhoff, A.K. and Moser, M. 2004. Identification and embryonic expression of a new AP-2 transcription factor, AP-2 ϵ . *Dev. Dyn.* 231: 128-135.
- Eckert, D., Buhl, S., Weber, S., Jäger, R. and Schorle, H. 2005. The AP-2 family of transcription factors. *Genome Biol.* 6: 246.
- Wenke, A.K., Rothhammer, T., Moser, M. and Bosserhoff, A.K. 2006. Regulation of integrin α 10 expression in chondrocytes by the transcription factors AP-2 ϵ and Ets-1. *Biochem. Biophys. Res. Commun.* 345: 495-501.
- Orso, F., Fassetta, M., Penna, E., Solero, A., De Filippo, K., Sismondì, P., De Bortoli, M. and Taverna, D. 2007. The AP-2 α transcription factor regulates tumor cell migration and apoptosis. *Adv. Exp. Med. Biol.* 604: 87-95.
- Orso, F., Penna, E., Cimino, D., Astanina, E., Maione, F., Valdembrì, D., Giraudo, E., Serini, G., Sismondì, P., De Bortoli, M. and Taverna, D. 2008. AP-2 α and AP-2 γ regulate tumor progression via specific genetic programs. *FASEB J.* 22: 2702-2714.

CHROMOSOMAL LOCATION

Genetic locus: TFAP2E (human) mapping to 1p34.3.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

AP-2 ϵ (T-25) is a purified rabbit polyclonal antibody raised against AP-2 ϵ of human origin.

PRODUCT

Each vial contains 100 μ g IgG in 1.0 ml PBS with < 0.1% sodium azide, 0.1% gelatin and < 0.02% sucrose.

APPLICATIONS

AP-2 ϵ (T-25) is recommended for detection of AP-2 ϵ of human and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Suitable for use as control antibody for AP-2 ϵ siRNA (h): sc-88389, AP-2 ϵ shRNA Plasmid (h): sc-88389-SH and AP-2 ϵ shRNA (h) Lentiviral Particles: sc-88389-V.

Molecular Weight of AP-2 ϵ : 46 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



AP-2 ϵ (T-25): sc-101887. Western blot analysis of AP-2 ϵ expression in transfected 293T whole cell lysate.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.