

FOXB1 (R-21): sc-101958

BACKGROUND

The Forkhead-box (FOX) genes comprise a superfamily of at least 43 members that encode proteins which are involved in transcriptional regulation and may be associated with the pathogenesis of various cancers. FOXB1 (forkhead box B1), also known as FKH5 or HFKH-5, and FOXB2 (forkhead box B2) are members of the FOX family and each contain one forkhead DNA-binding domain. Both FOXB1 and FOXB2 localize to the nucleus where they are thought to function as transcription factors that can bind to DNA via their forkhead domains. In mice, defects in the gene encoding FOXB1 are associated with retarded development of the central nervous system (CNS), suggesting that FOXB1 may play a role in CNS organization and function.

REFERENCES

- Weigel, D. and Jäckle, H. 1990. The fork head domain: a novel DNA binding motif of eukaryotic transcription factors? *Cell* 63: 455-456.
- Kaestner, K.H., Schütz, G. and Monaghan, A.P. 1996. Expression of the winged helix genes fkh-4 and fkh-5 defines domains in the central nervous system. *Mech. Dev.* 55: 221-230.
- Wehr, R., Mansouri, A., de Maeyer, T. and Gruss, P. 1997. Fkh5-deficient mice show dysgenesis in the caudal midbrain and hypothalamic mammillary body. *Development* 124: 4447-4456.
- Alvarez-Bolado, G., Cecconi, F., Wehr, R. and Gruss, P. 1999. The fork head transcription factor Fkh5/Mf3 is a developmental marker gene for superior colliculus layers and derivatives of the hindbrain somatic afferent zone. *Brain Res. Dev. Brain Res.* 112: 205-215.
- Alvarez-Bolado, G., Zhou, X., Cecconi, F. and Gruss, P. 2000. Expression of Foxb1 reveals two strategies for the formation of nuclei in the developing ventral diencephalon. *Dev. Neurosci.* 22: 197-206.
- Kloetzli, J.M., Fontaine-Glover, I.A., Brown, E.R., Kuo, M. and Labosky, P.A. 2001. The winged helix gene, Foxb1, controls development of mammary glands and regions of the CNS that regulate the milk-ejection reflex. *Genesis* 29: 60-71.
- Pohl, B.S., Knöchel, S., Dillinger, K. and Knöchel, W. 2002. Sequence and expression of FoxB2 (XFD-5) and Foxl1c (XFD-10) in *Xenopus* embryogenesis. *Mech. Dev.* 117: 283-287.
- Katoh, M. and Katoh, M. 2004. Human FOX gene family (Review). *Int. J. Oncol.* 25: 1495-1500.
- Radyushkin, K., Anokhin, K., Meyer, B.I., Jiang, Q., Alvarez-Bolado, G. and Gruss, P. 2005. Genetic ablation of the mammillary bodies in the Foxb1 mutant mouse leads to selective deficit of spatial working memory. *Eur. J. Neurosci.* 21: 219-229.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: FOXB1 (human) mapping to 15q22.2; Foxb1 (mouse) mapping to 9 D.

SOURCE

FOXB1 (R-21) is a purified rabbit polyclonal antibody raised against FOXB1 of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml PBS with < 0.1% sodium azide, 0.1% gelatin and < 0.02% sucrose.

APPLICATIONS

FOXB1 (R-21) is recommended for detection of FOXB1 of mouse, human, dog and zebrafish origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for FOXB1 siRNA (h): sc-90235, FOXB1 siRNA (m): sc-145219, FOXB1 shRNA Plasmid (h): sc-90235-SH, FOXB1 shRNA Plasmid (m): sc-145219-SH, FOXB1 shRNA (h) Lentiviral Particles: sc-90235-V and FOXB1 shRNA (m) Lentiviral Particles: sc-145219-V.


Molecular Weight of FOXB1: 35 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



FOXB1 (R-21): sc-101958. Western blot analysis of FOXB1 expression in Hep G2 whole cell lysate.



FOXB1 (R-21): sc-101958. Western blot analysis of FOXB1 expression in Hep G2 nuclear extract.

RESEARCH USE

For research use only, not for use in diagnostic procedures.