

RAD54L2 (G-18): sc-103135

BACKGROUND

Adrenergic receptors (ARs) include four general types (a1, a2, b1 and b2) that are found in different target tissues and differ in their affinities and responses to various agonists and antagonists. The coupling of ARs to specific intracellular effectors is mediated through diverse heterotrimeric G proteins. ARs play a critical role in the development of prostate cancer, and transcriptional activity of AR is partly regulated by coregulatory proteins. RAD54L2 (RAD54-like 2), also known as ARIP4 (androgen receptor-interacting protein 4), HSPC325 or SRISNF2L, is a 1,467 amino acid nuclear protein belonging to the SNF2/RAD54 helicase family that consists of one helicase ATP-binding domain and a helicase C-terminal domain. RAD54L2 is a DNA helicase that regulates androgen receptor (AR)-dependent transactivation in a promoter-dependent manner. RAD54L2 is posttranslationally sumoylated or phosphorylated upon DNA damage.

REFERENCES

- Jänne, O.A., et al. 2000. Androgen-receptor-interacting nuclear proteins. *Biochem. Soc. Trans.* 28: 401-405.
- Rouleau, N., et al. 2002. Novel ATPase of SNF2-like protein family interacts with androgen receptor and modulates androgen-dependent transcription. *Mol. Biol. Cell* 13: 2106-2119.
- Linja, M.J., et al. 2004. Expression of androgen receptor coregulators in prostate cancer. *Clin. Cancer Res.* 10: 1032-1040.
- Domanskyi, A., et al. 2006. Biochemical characterization of androgen receptor-interacting protein 4. *Biochem. J.* 393: 789-795.
- Rosendorff, A., et al. 2006. NXP-2 association with SUMO-2 depends on lysines required for transcriptional repression. *Proc. Natl. Acad. Sci. USA* 103: 5308-5313.
- Domanskyi, A., et al. 2007. Expression and localization of androgen receptor-interacting protein-4 in the testis. *Am. J. Physiol. Endocrinol. Metab.* 292: E513-E522.
- Zhang, F.P., et al. 2007. An adenosine triphosphatase of the sucrose non-fermenting 2 family, androgen receptor-interacting protein 4, is essential for mouse embryonic development and cell proliferation. *Mol. Endocrinol.* 21: 1430-1442.
- Urbanucci, A., et al. 2008. Androgen regulation of the androgen receptor coregulators. *BMC Cancer* 8: 219.
- Ogawa, H., et al. 2009. Transcriptional suppression by transient recruitment of ARIP4 to sumoylated nuclear receptor Ad4BP/SF-1. *Mol. Biol. Cell* 20: 4235-4245.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: RAD54L2 (human) mapping to 3p21.2.

SOURCE

RAD54L2 (G-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of RAD54L2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-103135 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

RAD54L2 (G-18) is recommended for detection of RAD54L2 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other Rad family members.


RAD54L2 (G-18) is also recommended for detection of RAD54L2 in additional species, including equine, canine and avian.

Suitable for use as control antibody for RAD54L2 siRNA (h): sc-78533, RAD54L2 shRNA Plasmid (h): sc-78533-SH and RAD54L2 shRNA (h) Lentiviral Particles: sc-78533-V.

Molecular Weight of RAD54L2: 180 kDa.

Positive Controls: RAD54L2 (h): 293T Lysate: sc-175676.

DATA



RAD54L2 (G-18): sc-103135. Western blot analysis of RAD54L2 expression in non-transfected: sc-117752 (A) and human RAD54L2 transfected: sc-175676 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.