



DCAMKL2 (S-13): sc-103449

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. A group of proteins that are intimately involved in this process are the serine/threonine (Ser/Thr) protein kinases. DCAMKL2, also known as DCLK2 (Doublecortin-like kinase 2), DCK2, DCDC3 or DCDC3B, is a 766 amino acid protein that contains one protein kinase domain and two Doublecortin domains. One of several members of the Ser/Thr protein kinase family, DCAMKL2 functions to catalyze the ATP-dependent attachment of a phospho residue to target proteins, an event that may play a role in various cellular processes. Multiple isoforms of DCAMKL2 exist due to alternative splicing events.

REFERENCES

- Bairoch, A. and Claverie, J.M. 1988. Sequence patterns in protein kinases. *Nature*. 331: 22.
- Hanks, S.K., Quinn, A.M. and Hunter, T. 1988. The protein kinase family: conserved features and deduced phylogeny of the catalytic domains. *Science* 241: 42-52.
- Hanks, S.K. and Quinn, A.M. 1991. Protein kinase catalytic domain sequence database: identification of conserved features of primary structure and classification of family members. *Meth. Enzymol.* 200: 38-62.
- Sapir, T., Horesh, D., Caspi, M., Atlas, R., Burgess, H.A., Wolf, S.G., Francis, F., Chelly, J., Elbaum, M., Pietrokovski, S. and Reiner, O. 2000. Doublecortin mutations cluster in evolutionarily conserved functional domains. *Hum. Mol. Genet.* 9: 703-712.
- Edelman, A.M., Kim, W.Y., Higgins, D., Goldstein, E.G., Oberdoerster, M. and Sigurdson, W. 2005. Doublecortin kinase-2, a novel Doublecortin-related protein kinase associated with terminal segments of axons and dendrites. *J. Biol. Chem.* 280: 8531-8543.
- Ohmae, S., Takemoto-Kimura, S., Okamura, M., Adachi-Morishima, A., Nonaka, M., Fuse, T., Kida, S., Tanji, M., Furuyashiki, T., Arakawa, Y., Narumiya, S., Okuno, H. and Bito, H. 2006. Molecular identification and characterization of a family of kinases with homology to Ca²⁺/calmodulin-dependent protein kinases I/IV. *J. Biol. Chem.* 281: 20427-20439.
- Koizumi, H. 2007. Physiological function of the DCX family in brain development. *Seikagaku* 79: 1134-1139.
- Tuy, F.P., Saillour, Y., Kappeler, C., Chelly, J. and Francis, F. 2008. Alternative transcripts of DCLK1 and DCLK2 and their expression in Doublecortin knock-out mice. *Dev. Neurosci.* 30: 171-186.

CHROMOSOMAL LOCATION

Genetic locus: DCLK2 (human) mapping to 4q31.3; Dclk2 (mouse) mapping to 3 F1.

SOURCE

DCAMKL2 (S-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of DCAMKL2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-103449 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

DCAMKL2 (S-13) is recommended for detection of DCAMKL2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with DCAMKL1 or DCAMKL3.

Suitable for use as control antibody for DCAMKL2 siRNA (h): sc-105272, DCAMKL2 siRNA (m): sc-105273, DCAMKL2 shRNA Plasmid (h): sc-105272-SH, DCAMKL2 shRNA Plasmid (m): sc-105273-SH, DCAMKL2 shRNA (h) Lentiviral Particles: sc-105272-V and DCAMKL2 shRNA (m) Lentiviral Particles: sc-105273-V.

Molecular Weight of DCAMKL2: 84 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.