

# ACOT7 (S-15): sc-104029

## BACKGROUND

Acyl-CoA thioesterases (ACOTs) are a group of enzymes that catalyze the hydrolysis of acyl-CoA to form coenzyme A (CoA) and a free fatty acid. Through their catalytic activity, ACOTs are able to regulate the level of fatty acids and acyl-CoAs within the cell. ACOT7 (acyl-CoA thioesterase 7), also known as BACH (brain acyl-CoA hydrolase), LACH or CTE-II, is a 380 amino acid protein that is expressed as 6 alternatively spliced isoforms which localize to either the cytoplasm or the mitochondria. Functioning as a homodimer that contains 2 acyl coenzyme A hydrolase domains, ACOT7 plays an important role in regulating acyl-CoA levels within the body and is thought to specifically participate in proper brain physiology and function. Decreased ACOT7 expression may be associated with mesial temporal lobe epilepsy, a form of focal epilepsy that is characterized by simple or complex seizures.

## REFERENCES

1. Yamada, J., Kurata, A., Hirata, M., Taniguchi, T., Takama, H., Furihata, T., Shiratori, K., Iida, N., Takagi-Sakuma, M., Watanabe, T., Kurosaki, K., Endo, T. and Suga, T. 1999. Purification, molecular cloning, and genomic organization of human brain long-chain acyl-CoA hydrolase. *J. Biochem.* 126: 1013-1019.
2. Yamada, J., Kuramochi, Y., Takagi, M., Watanabe, T. and Suga, T. 2002. Human brain acyl-CoA hydrolase isoforms encoded by a single gene. *Biochem. Biophys. Res. Commun.* 299: 49-56.
3. Kuramochi, Y., Takagi-Sakuma, M., Kitahara, M., Emori, R., Asaba, Y., Sakaguchi, R., Watanabe, T., Kuroda, J., Hiratsuka, K., Nagae, Y., Suga, T. and Yamada, J. 2002. Characterization of mouse homolog of brain acyl-CoA hydrolase: molecular cloning and neuronal localization. *Brain Res. Mol. Brain Res.* 98: 81-92.
4. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602587. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Yang, J.W., Czech, T., Yamada, J., Cszasz, E., Baumgartner, C., Slavic, I. and Lubec, G. 2004. Aberrant cytosolic acyl-CoA thioester hydrolase in hippocampus of patients with mesial temporal lobe epilepsy. *Amino Acids* 27: 269-275.
6. Yamada, J. 2005. Long-chain acyl-CoA hydrolase in the brain. *Amino Acids* 28: 273-278.
7. Hunt, M.C., Yamada, J., Maltais, L.J., Wright, M.W., Podesta, E.J. and Alexson, S.E. 2005. A revised nomenclature for mammalian acyl-CoA thioesterases/hydrolases. *J. Lipid Res.* 46: 2029-2032.
8. Serek, R., Forwood, J.K., Hume, D.A., Martin, J.L. and Kobe, B. 2006. Crystallization of the C-terminal domain of the mouse brain cytosolic long-chain acyl-CoA thioesterase. *Acta Crystallogr. Sect. F Struct. Biol. Cryst. Commun.* 62: 133-135.

## CHROMOSOMAL LOCATION

Genetic locus: ACOT7 (human) mapping to 1p36.31; Acot7 (mouse) mapping to 4 E2.

## SOURCE

ACOT7 (S-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of ACOT7 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-104029 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

ACOT7 (S-15) is recommended for detection of ACOT7 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

ACOT7 (S-15) is also recommended for detection of ACOT7 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for ACOT7 siRNA (h): sc-88501, ACOT7 siRNA (m): sc-105035, ACOT7 shRNA Plasmid (h): sc-88501-SH, ACOT7 shRNA Plasmid (m): sc-105035-SH, ACOT7 shRNA (h) Lentiviral Particles: sc-88501-V and ACOT7 shRNA (m) Lentiviral Particles: sc-105035-V.

Molecular Weight of ACOT7 isoforms B/A-X/A-Xi: 42/27/31 kDa.

Molecular Weight of ACOT7 isoforms A/C/D: 37/39/37 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **ACOT7 (C-2): sc-376808** or **ACOT7 (B-4): sc-376692**, our highly recommended monoclonal alternatives to ACOT7 (S-15).