

Arp-T2 (K-12): sc-104073

BACKGROUND

Actin-related proteins are classified into Arp subclasses according to their amino acid sequence similarity to Actin. Both Arps and Actin proteins have an ATPase domain, which catalyzes the decomposition of adenosine triphosphate (ATP) into adenosine diphosphate (ADP) and a free phosphate ion to release energy. Arp-T2 (Actin-related protein T2), also known as Arp-M2, is a 377-amino acid structural cytoskeleton protein. Like other Arp family members, Arp-T2 contains an Actin-like ATP/ADP-binding pocket and two nuclear export signals. Although Arp-T2 shares 43% sequence identity with β -Actin, it contains several cysteine residues that are not found in other Actin family members. With significant expression in testis, Arp-T2 may play a role in spermatid formation.

REFERENCES

1. Longo, F.J., Krohne, G. and Franke, W.W. 1987. Basic proteins of the perinuclear theca of mammalian spermatozoa and spermatids: a novel class of cytoskeletal elements. *J. Cell Biol.* 105: 1105-1120.
2. Lecuyer, C., Dacheux, J.L., Hermand, E., Mazeman, E., Rousseaux, J. and Rousseaux-Prevost, R. 2000. Actin-binding properties and co-localization with Actin during spermiogenesis of mammalian sperm calicin. *Biol. Reprod.* 63: 1801-1810.
3. Harata, M., Nishimori, K. and Hatta, S. 2001. Identification of two cDNAs for human Actin-related proteins (Arps) that have remarkable similarity to conventional Actin. *Biochim. Biophys. Acta* 1522: 130-133.
4. Heid, H., Figge, U., Winter, S., Kuhn, C., Zimbelmann, R. and Franke, W. 2002. Novel Actin-related proteins Arp-T1 and Arp-T2 as components of the cytoskeletal calyx of the mammalian sperm head. *Exp. Cell Res.* 279: 177-187.
5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 608535. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Rousseaux-Prevost, R., Lecuyer, C., Drobecq, H., Sergheraert, C., Dacheux, J.L. and Rousseaux, J. 2003. Characterization of boar sperm cytoskeletal cylicin II as an Actin-binding protein. *Biochem. Biophys. Res. Commun.* 303: 182-189.
7. Hara, Y., Yamagata, K., Oguchi, K. and Baba, T. 2008. Nuclear localization of profilin III-ArpM1 complex in mouse spermiogenesis. *FEBS Lett.* 582: 2998-3004.

CHROMOSOMAL LOCATION

Genetic locus: ACTRT2 (human) mapping to 1p36.32.

SOURCE

Arp-T2 (K-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Arp-T2 of human origin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-104073 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Arp-T2 (K-12) is recommended for detection of Arp-T2 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Arp-T2 siRNA (h): sc-88470, Arp-T2 shRNA Plasmid (h): sc-88470-SH and Arp-T2 shRNA (h) Lentiviral Particles: sc-88470-V.

Molecular Weight (predicted) of Arp-T2: 42 kDa.


Molecular Weight (observed) of Arp-T2: 37/43 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227, T98G cell lysate: sc-2294 or RT-4 whole cell lysate: sc-364257.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Arp-T2 (K-12): sc-104073. Western blot analysis of Arp-T2 expression in T98G (A), Hep G2 (B) and RT-4 (C) whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.