

GATAD1 (C-16): sc-104262

BACKGROUND

GATAD1 (GATA zinc finger domain-containing protein 1), also known as ODAG (ocular development-associated gene protein), is a 269 amino acid protein that is involved in early ocular development. Expressed highly in postnatal eye tissue, GATAD1 is associated with formation of the lens and its surrounding structures, suggesting a possible role in the transformation of ocular tissues into a working eye. GATAD1 expression declines dramatically after the early stages of development. GATAD1 contains 1 GATA-type zinc finger which functions as a DNA-binding domain. Additionally, GATAD1 gene expression is amplified in certain cancerous cells, suggesting that it may be involved in carcinogenesis.

REFERENCES

1. Tsuruga, T., Kanamoto, T., Kato, T., Yamashita, H., Miyagawa, K. and Mishima, H.K. 2002. Ocular development-associated gene (ODAG), a novel gene highly expressed in ocular development. *Gene* 290: 125-130.
2. Pang, E.Y., Bai, A.H., To, K.F., Sy, S.M., Wong, N.L., Lai, P.B., Squire, J.A. and Wong, N. 2007. Identification of PFTAIRE protein kinase 1, a novel cell division cycle-2 related gene, in the motile phenotype of hepatocellular carcinoma cells. *Hepatology* 46: 436-445.
3. Agelopoulos, K., Kersting, C., Korsching, E., Schmidt, H., Kuijper, A., August, C., Wülfing, P., Tio, J., Boecker, W., van Diest, P.J., Brandt, B. and Buerger, H. 2007. EGFR amplification specific gene expression in phyllodes tumours of the breast. *Cell. Oncol.* 29: 443-451.
4. Nagel, S., Leich, E., Quentmeier, H., Meyer, C., Kaufmann, M., Drexler, H.G., Zettl, A., Rosenwald, A. and Macleod, R.A. 2008. Amplification at 7q22 targets cyclin-dependent kinase 6 in T cell lymphoma. *Leukemia* 22: 387-392.
5. Sasaki, T., Watanabe, W., Muranishi, Y., Kanamoto, T., Aihara, M., Miyazaki, K., Tamura, H., Saeki, T., Oda, H., Souchelnytskyi, N., Souchelnytskyi, S., Aoyama, H., Honda, Z.I., Furukawa, T., Mishima, H.K., Kiuchi, Y. and Honda, H. 2008. Overexpression of ocular development-associated gene (ODAG) in the eye induced elevated intraocular pressure, optic nerve atrophy, and impaired retinal development. *Invest. Ophthalmol. Vis. Sci.* 50: 242-248.

CHROMOSOMAL LOCATION

Genetic locus: GATAD1 (human) mapping to 7q21.2; Gatad1 (mouse) mapping to 5 A1.

SOURCE

GATAD1 (C-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of GATAD1 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-104262 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

GATAD1 (C-16) is recommended for detection of GATAD1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with family members GATAD2A or GATAD2B.

Suitable for use as control antibody for GATAD1 siRNA (h): sc-89708, GATAD1 siRNA (m): sc-145341, GATAD1 shRNA Plasmid (h): sc-89708-SH, GATAD1 shRNA Plasmid (m): sc-145341-SH, GATAD1 shRNA (h) Lentiviral Particles: sc-89708-V and GATAD1 shRNA (m) Lentiviral Particles: sc-145341-V.

Molecular Weight of GATAD1: 29 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200 or NIH/3T3 whole cell lysate: sc-2210.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.