

HS3ST1 (C-15): sc-104313

BACKGROUND

Heparan sulfate structures, which are responsible for executing multiple biologic activities, are generated and regulated by heparan sulfate biosynthetic enzymes. HS3ST1 (heparan sulfate (glucosamine) 3-O-sulfotransferase 1), also known as 3OST or 3OST1, is a 307 amino acid protein that localizes to the lumen of the Golgi apparatus and belongs to the heparan sulfate biosynthetic enzyme family. Expressed at high levels in kidney and brain, and present at lower levels in lung, heart and placenta, HS3ST1 exhibits both anticoagulant heparan sulfate conversion activity and heparan sulfate glucosaminyl 3-O-sulfotransferase activity and specifically catalyzes the rate-limiting step in the synthesis of anticoagulant heparan. The gene encoding HS3ST1 maps to human chromosome 4, which encodes nearly 6% of the human genome and has the largest gene deserts (regions of the genome with no protein encoding genes) of all of the human chromosomes.

REFERENCES

1. Razi, N. and Lindahl, U. 1995. Biosynthesis of heparin/heparan sulfate. The D-glucosaminyl 3-O-sulfotransferase reaction: target and inhibitor saccharides. *J. Biol. Chem.* 270: 11267-11275.
2. Liu, J., et al. 1996. Purification of heparan sulfate D-glucosaminyl 3-O-sulfotransferase. *J. Biol. Chem.* 271: 27072-27082.
3. Shworak, N.W., et al. 1997. Molecular cloning and expression of mouse and human cDNAs encoding heparan sulfate D-glucosaminyl 3-O-sulfotransferase. *J. Biol. Chem.* 272: 28008-28019.
4. Shworak, N.W., et al. 1999. Multiple isoforms of heparan sulfate D-glucosaminyl 3-O-sulfotransferase. Isolation, characterization, and expression of human cDNAs and identification of distinct genomic loci. *J. Biol. Chem.* 274: 5170-5184.

CHROMOSOMAL LOCATION

Genetic locus: HS3ST1 (human) mapping to 4p15.33; Hs3st1 (mouse) mapping to 5 B3.

SOURCE

HS3ST1 (C-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of HS3ST1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-104313 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

HS3ST1 (C-15) is recommended for detection of HS3ST1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other HS3ST family members.

HS3ST1 (C-15) is also recommended for detection of HS3ST1 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for HS3ST1 siRNA (h): sc-88999, HS3ST1 siRNA (m): sc-105541, HS3ST1 shRNA Plasmid (h): sc-88999-SH, HS3ST1 shRNA Plasmid (m): sc-105541-SH, HS3ST1 shRNA (h) Lentiviral Particles: sc-88999-V and HS3ST1 shRNA (m) Lentiviral Particles: sc-105541-V.


Molecular Weight of HS3ST1: 46 kDa.

Positive Controls: mouse kidney extract: sc-2255.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941. 4) Immunohistochemistry: use ImmunoCruz™: sc-2053 or ABC: sc-2023 goat IgG Staining Systems.

DATA



HS3ST1 (C-15): sc-104313. Western blot analysis of HS3ST1 expression in mouse kidney tissue extract.



HS3ST1 (C-15): sc-104313. Immunoperoxidase staining of formalin fixed, paraffin-embedded human hippocampus tissue showing cytoplasmic and nuclear staining of glial cells and cytoplasmic staining of neuronal cells.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.