

# MRP-S11 (P-15): sc-104407

## BACKGROUND

Mitochondrial ribosomes consist of a large 39S subunit and a small 28S subunit, both of which are comprised of multiple mitochondrial ribosomal proteins (MRPs) that are encoded by nuclear genes and are essential for protein synthesis within mitochondria. MRP-S11 (mitochondrial ribosomal protein S11) also known as Cervical cancer proto-oncogene 2 protein, is a 194 amino acid protein that localizes to the mitochondrion, where it exists as a component of the 28S ribosomal subunit and works in conjunction with other MRPs to mediate protein synthesis. Mutations in the MRP-S16 and MRP-S22 genes also cause a loss of MRP-S11 expression. The phenotype of these gene mutations is severe fatal respiratory chain dysfunction due to impaired translation of mitochondrial mRNAs. There are three isoforms of MRP-S11 that are produced as a result of alternative splicing events.

## REFERENCES

1. Kenmochi, N., Suzuki, T., Uechi, T., Magoori, M., Kuniba, M., Higa, S., Watanabe, K. and Tanaka, T. 2001. The human mitochondrial ribosomal protein genes: mapping of 54 genes to the chromosomes and implications for human disorders. *Genomics* 77: 65-70.
2. Cavdar Koc, E., Burkhart, W., Blackburn, K., Moseley, A. and Spremulli, L.L. 2001. The small subunit of the mammalian mitochondrial ribosome. Identification of the full complement of ribosomal proteins present. *J. Biol. Chem.* 276: 19363-19374.
3. Suzuki, T., Terasaki, M., Takemoto-Hori, C., Hanada, T., Ueda, T., Wada, A. and Watanabe, K. 2001. Proteomic analysis of the mammalian mitochondrial ribosome. Identification of protein components in the 28 S small subunit. *J. Biol. Chem.* 276: 33181-33195.
4. Zhang, Z. and Gerstein, M. 2003. Identification and characterization of over 100 mitochondrial ribosomal protein pseudogenes in the human genome. *Genomics* 81: 468-480.
5. O'Brien, T.W., O'Brien, B.J. and Norman, R.A. 2005. Nuclear MRP genes and mitochondrial disease. *Gene* 354: 147-151.
6. Smits, P., Smeitink, J.A., van den Heuvel, L.P., Huynen, M.A. and Ettema, T.J. 2007. Reconstructing the evolution of the mitochondrial ribosomal proteome. *Nucleic Acids Res.* 35: 4686-4703.
7. Cavaliere, D., Dolara, P., Mini, E., Luceri, C., Castagnini, C., Toti, S., Maciag, K., De Filippo, C., Nobili, S., Morganti, M., Napoli, C., Tonini, G., Baccini, M., Biggeri, A., Tonelli, F., Valanzano, R., Orlando, C., Gelmini, S., Cianchi, F., Messerini, L. and Luzzatto, L. 2007. Analysis of gene expression profiles reveals novel correlations with the clinical course of colorectal cancer. *Oncol. Res.* 16: 535-548.
8. Emdadul Haque, M., Grasso, D., Miller, C., Spremulli, L.L. and Saada, A. 2008. The effect of mutated mitochondrial ribosomal proteins S16 and S22 on the assembly of the small and large ribosomal subunits in human mitochondria. *Mitochondrion* 8: 254-261.
9. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 611977. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

## CHROMOSOMAL LOCATION

Genetic locus: MRPS11 (human) mapping to 15q25.3; Mrps11 (mouse) mapping to 7 D3.

## SOURCE

MRP-S11 (P-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of MRP-S11 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-104407 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

MRP-S11 (P-15) is recommended for detection of MRP-S11 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with MRP-S11-3.

Suitable for use as control antibody for MRP-S11 siRNA (h): sc-90071, MRP-S11 siRNA (m): sc-106250, MRP-S11 shRNA Plasmid (h): sc-90071-SH, MRP-S11 shRNA Plasmid (m): sc-106250-SH, MRP-S11 shRNA (h) Lentiviral Particles: sc-90071-V and MRP-S11 shRNA (m) Lentiviral Particles: sc-106250-V.

Molecular Weight of MRP-S11 isoforms: 21/17 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker<sup>™</sup> compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz<sup>™</sup> Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.