

# FXR1 (C-18): sc-10554

## BACKGROUND

Fragile X syndrome is the most frequent form of inherited mental retardation and is the result of transcriptional silencing of the FMR1 gene on the X chromosome. The FMR1 gene contains a distinct CpG dinucleotide repeat located in the 5' untranslated region of the gene, and in the fragile X syndrome this tandem repeat is substantially amplified, and subjected to extensive methylation and enhanced transcriptional silencing. The FMR1 protein (or FMRP) is an RNA-binding protein that associates with polyribosomes and is a likely component of a messenger ribonuclear protein (mRNP) particle. FMR1 contains several features that are characteristics of RNA-binding proteins, including two hnRNPK homology (KH) domains and an RGG amino acid motif (RGG box). FMR1 can also interact with two fragile X syndrome related factors, FXR1 and FXR2, and these proteins form heterodimers through their N-terminal coil-coiled domains. FMR1 localizes to both the nucleus and the cytoplasm, and since it contains both a nuclear localization signal and a nuclear export signal it is also implicated in the nucleocytoplasmic transport of mRNAs.

## CHROMOSOMAL LOCATION

Genetic locus: FXR1 (human) mapping to 3q26.33; Fxr1 (mouse) mapping to 3 A3.

## SOURCE

FXR1 (C-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of FXR1 of human origin.

## PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-10554 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## APPLICATIONS

FXR1 (C-18) is recommended for detection of FXR1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

FXR1 (C-18) is also recommended for detection of FXR1 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for FXR1 siRNA (h): sc-35423, FXR1 siRNA (m): sc-35424, FXR1 shRNA Plasmid (h): sc-35423-SH, FXR1 shRNA Plasmid (m): sc-35424-SH, FXR1 shRNA (h) Lentiviral Particles: sc-35423-V and FXR1 shRNA (m) Lentiviral Particles: sc-35424-V.


Molecular Weight of FXR1: 78 kDa.

Positive Controls: F9 cell lysate: sc-2245 or MCF7 whole cell lysate: sc-2206.


## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



FXR1 (C-18): sc-10554. Western blot analysis of FXR1 expression in F9 (A) and MCF7 (B) whole cell lysate.



FXR1 (C-18): sc-10554. Immunofluorescence staining of methanol-fixed F9 cells showing cytoplasmic localization.

## SELECT PRODUCT CITATIONS

- Otsuka, Y., et al. 2009. Identification of a cytoplasmic complex that adds a cap onto 5'-monophosphate RNA. *Mol. Cell. Biol.* 29: 2155-2167.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.

**MONOS**  
Satisfaction  
Guaranteed

Try **FXR1 (B-2): sc-374148** or **FXR1 (C-4): sc-515398**, our highly recommended monoclonal alternatives to FXR1 (C-18). Also, for AC, HRP, FITC, PE, Alexa Fluor® 488 and Alexa Fluor® 647 conjugates, see **FXR1 (B-2): sc-374148**.