

FXR2 (D-18): sc-10556

BACKGROUND

Fragile X syndrome is the most frequent form of inherited mental retardation and is the result of transcriptional silencing of the FMR1 gene on the X chromosome. The FMR1 gene contains a distinct CpG dinucleotide repeat located in the 5' untranslated region of the gene, and in the fragile X syndrome this tandem repeat is substantially amplified, and subjected to extensive methylation and enhanced transcriptional silencing. The FMR1 protein (or FMRP) is an RNA-binding protein that associates with polyribosomes and is a likely component of a messenger ribonuclear protein (mRNP) particle. FMR1 contains several features that are characteristics of RNA-binding proteins, including two hnRNPK homology (KH) domains and an RGG amino acid motif (RGG box). FMR1 can also interact with two fragile X syndrome related factors, FXR1 and FXR2, and these proteins form heterodimers through their N-terminal coil-coiled domains. FMR1 localizes to both the nucleus and the cytoplasm, and since it contains both a nuclear localization signal and a nuclear export signal it is also implicated in the nucleo-cytoplasmic transport of mRNAs.

REFERENCES

1. Verkerk, A.J., et al. 1991. Identification of a gene (FMR1) containing a CGG repeat coincident with a breakpoint cluster region exhibiting length variation in fragile X syndrome. *Cell* 65: 905-914.
2. Pieretti, M., et al. 1991. Absence of expression of the FMR1 gene in fragile X syndrome. *Cell* 66: 817-22.
3. Matunis, M.J., et al. 1992. Characterization and primary structure of the poly(C)-binding heterogeneous nuclear ribo-nucleoprotein complex K protein. *Mol. Cell. Biol.* 12: 164-171.
4. De Boulle, K., et al. 1993. A point mutation in the FMR1 gene associated with fragile X mental retardation. *Nat. Genet.* 3: 31-35.
5. Zhang, Y., et al. 1995. The fragile X mental retardation syndrome protein interacts with novel homologs FXR1 and FXR2. *EMBO J.* 14: 5358-5366.
6. Eberhart, D.E., et al. 1996. The fragile X mental retardation protein is a ribonucleoprotein containing both nuclear localization and nuclear export signals. *Hum. Mol. Genet.* 5: 1083-1091.
7. Ceman, S., et al. 1999. Isolation of an FMRP-associated messenger ribonucleoprotein particle and identification of nucleolin and the fragile X related proteins as components of the complex. *Mol. Cell. Biol.* 19: 7925-7932.

CHROMOSOMAL LOCATION

Genetic locus: FXR2 (human) mapping to 17p13.1; Fxr2h (mouse) mapping to 11 B3.

SOURCE

FXR2 (D-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of FXR2 of mouse origin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-10556 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

FXR2 (D-18) is recommended for detection of FXR2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


FXR2 (D-18) is also recommended for detection of FXR2 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for FXR2 siRNA (h): sc-37895, FXR2 siRNA (m): sc-37896, FXR2 shRNA Plasmid (h): sc-37895-SH, FXR2 shRNA Plasmid (m): sc-37896-SH, FXR2 shRNA (h) Lentiviral Particles: sc-37895-V and FXR2 shRNA (m) Lentiviral Particles: sc-37896-V.

Molecular Weight of FXR2: 90 kDa.

Positive Controls: FXR2 (h3): 293T Lysate: sc-117244, IMR-32 cell lysate: sc-2409 or F9 cell lysate: sc-2245.

DATA



FXR2 (D-18): sc-10556. Western blot analysis of FXR2 expression in non-transfected: sc-117752 (A) and human FXR2 transfected: sc-117244 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **FXR2 (1G2): sc-32266** or **FXR2 (G-2): sc-376963**, our highly recommended monoclonal alternatives to FXR2 (D-18).