

CASKIN1 (C-13): sc-107180

BACKGROUND

CASK is a multi-domain scaffolding protein that is thought to be involved in signaling at neuronal synapses. CASKIN1 (CASK-interacting protein 1) is a 1,431 amino acid cytoplasmic protein that links CASK to downstream intracellular effectors. CASKIN1 contains six ANK repeats, two SAM (sterile α motif) domains and one SH3 domain. By binding the CaM kinase domain of CASK, CASKIN1 forms a ternary complex with CASK and a Veli protein, either 1, 2 or 3. This tripartite binds to the C-terminus of neurexin I, a neuronal cell surface protein involved in cell recognition and adhesion. X11 α , a CNS-specific protein that may be involved in the processing of the β -amyloid precursor protein, also binds the CASK and LIN7 complex and therefore competes with CASKIN1 for CASK binding.

REFERENCES

1. Tabuchi, K., Biederer, T., Butz, S. and Sudhof, T.C. 2002. CASK participates in alternative tripartite complexes in which Mint 1 competes for binding with CASKIN1, a novel CASK-binding protein. *J. Neurosci.* 22: 4264-4273.
2. Wang, G.S., Hong, C.J., Yen, T.Y., Huang, H.Y., Ou, Y., Huang, T.N., Jung, W.G., Kuo, T.Y., Sheng, M., Wang, T.F. and Hsueh, Y.P. 2004. Transcriptional modification by a CASK-interacting nucleosome assembly protein. *Neuron* 42: 113-128.
3. Boehm, D., Bacher, J. and Neumann, H.P. 2007. Gross genomic rearrangement involving the TSC2-PKD1 contiguous deletion syndrome: characterization of the deletion event by quantitative polymerase chain reaction deletion assay. *Am. J. Kidney Dis.* 49: e11-e21.
4. Online Mendelian Inheritance in Man, OMIM™. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 612184. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Bailey, K.A. and Aldinger, K.A. 2009. An X-linked microcephaly syndrome caused by disruptions of CASK implicates the CASK-TBR1-RELN pathway in human brain development. *Clin. Genet.* 75: 424-425.
6. Balázs, A., Csizmok, V., Buday, L., Rakács, M., Kiss, R., Bokor, M., Udupa, R., Tompa, K. and Tompa, P. 2009. High levels of structural disorder in scaffold proteins as exemplified by a novel neuronal protein, CASK-interactive protein1. *FEBS J.* 276: 3744-3756.
7. Wang, Z., Wang, B., Yang, L., Guo, Q., Aithmitti, N., Songyang, Z. and Zheng, H. 2009. Presynaptic and postsynaptic interaction of the amyloid precursor protein promotes peripheral and central synaptogenesis. *J. Neurosci.* 29: 10788-10801.
8. Sun, M., Liu, L., Zeng, X., Xu, M., Liu, L., Fang, M. and Xie, W. 2009. Genetic interaction between Neurexin and CAKI/CMG is important for synaptic function in *Drosophila* neuromuscular junction. *Neurosci. Res.* 64: 362-371.

CHROMOSOMAL LOCATION

Genetic locus: CASKIN1 (human) mapping to 16p13.3; Caskin1 (mouse) mapping to 17 A3.3.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

CASKIN1 (C-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of CASKIN1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-107180 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

CASKIN1 (C-13) is recommended for detection of CASKIN1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with family member CASKIN2.

CASKIN1 (C-13) is also recommended for detection of CASKIN1 in additional species, including bovine.

Suitable for use as control antibody for CASKIN1 siRNA (h): sc-93195, CASKIN1 siRNA (m): sc-142020, CASKIN1 shRNA Plasmid (h): sc-93195-SH, CASKIN1 shRNA Plasmid (m): sc-142020-SH, CASKIN1 shRNA (h) Lentiviral Particles: sc-93195-V and CASKIN1 shRNA (m) Lentiviral Particles: sc-142020-V.

Molecular Weight of CASKIN1: 180 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.