

TLR4 (H-80): sc-10741

BACKGROUND

Six human homologs of the *Drosophila* Toll receptor were initially identified based on their sequence similarities and designated Toll-like receptors (TLR). Toll receptors are involved in mediating dorsoventral polarization in the developing *Drosophila* embryo and also participate in the host immunity. The TLR family of proteins are characterized by a highly conserved Toll homology (TH) domain, which is essential for Toll-induced signal transduction. TLR1, as well as the other TLR family members, are type I transmembrane receptors that characteristically contain an extracellular domain consisting of several leucine-rich regions along with a single cytoplasmic Toll/IL-1R-like domain. TLR2 and TLR4 are activated in response to lipopolysaccharide (LPS) stimulation, which results in the activation and translocation of NF κ B and suggests that these receptors are involved in mediating inflammatory responses. Expression of TLR receptors is highest in peripheral blood leukocytes, macrophages, and monocytes. TLR6 is highly homologous to TLR1, sharing greater than 65% sequence identity, and, like other members of TLR family, it induces NF κ B signaling upon activation.

CHROMOSOMAL LOCATION

Genetic locus: TLR4 (human) mapping to 9q33.1; Tlr4 (mouse) mapping to 4 C1.

SOURCE

TLR4 (H-80) is a rabbit polyclonal antibody raised against amino acids 242-321 of TLR4 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

TLR4 (H-80) is recommended for detection of TLR4 of human and, to a lesser extent, mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for TLR4 siRNA (h): sc-40260, TLR4 siRNA (m): sc-40261, TLR4 siRNA (r): sc-156001, TLR4 shRNA Plasmid (h): sc-40260-SH, TLR4 shRNA Plasmid (m): sc-40261-SH, TLR4 shRNA Plasmid (r): sc-156001-SH, TLR4 shRNA (h) Lentiviral Particles: sc-40260-V, TLR4 shRNA (m) Lentiviral Particles: sc-40261-V and TLR4 shRNA (r) Lentiviral Particles: sc-156001-V.

Molecular Weight of TLR4: 95 kDa.

Molecular Weight of glycosylated TLR4: 120 kDa.

Positive Controls: HEL 92.1.7 cell lysate: sc-2270, THP-1 cell lysate: sc-2238 or HL-60 whole cell lysate: sc-2209.


STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.


RESEARCH USE

For research use only, not for use in diagnostic procedures.

DATA



TLR4 (H-80): sc-10741. Western blot analysis of TLR4 expression in HEL 92.1.7 (A) and HL-60 (B) whole cell lysates.



TLR4 (H-80): sc-10741. Western blot analysis of TLR4 expression in THP-1 whole cell lysate.

SELECT PRODUCT CITATIONS

- Medvedev, A.E., et al. 2002. Dysregulation of LPS-induced Toll-like receptor 4-MyD88 complex formation and IL-1 receptor-associated kinase 1 activation in endotoxin-tolerant cells. *J. Immunol.* 169: 5209-5216.
- Sienerth, A.R., et al. 2011. Polycomb group protein Bmi1 negatively regulates IL-10 expression in activated macrophages. *Immunol. Cell Biol.* 89: 812-816.
- Silva, M.A., et al. 2012. Increased bacterial translocation in gluten-sensitive mice is independent of small intestinal paracellular permeability defect. *Dig. Dis. Sci.* 57: 38-47.
- Pace, E., et al. 2012. β defensin-2 is reduced in central but not in distal airways of smoker COPD patients. *PLoS One* 7: e33601.
- Eiró, N., et al. 2012. Duodenal expression of toll-like receptors and interleukins are increased in both children and adult celiac patients. *Dig. Dis. Sci.* 57: 2278-2285.
- Pal, D., et al. 2012. Fetuin-A acts as an endogenous ligand of TLR4 to promote lipid-induced Insulin resistance. *Nat. Med.* 18: 1279-1285.
- Ismail, Y., et al. 2013. The effects of oral and enteric *Campylobacter concisus* strains on expression of TLR4, MD-2, TLR2, TLR5 and COX-2 in HT-29 cells. *PLoS ONE* 8: e56888.
- Ioannidis, I., et al. 2013. Toll-like receptor expression and induction of type I and type III interferons in primary airway epithelial cells. *J. Virol.* 87: 3261-3270.
- Mietto, B.S., et al. 2013. Lack of galectin-3 speeds Wallerian degeneration by altering TLR and pro-inflammatory cytokine expressions in injured sciatic nerve. *Eur. J. Neurosci.* 37: 1682-1690.
- Yao, L., et al. 2013. Toll-like receptor 4 mediates microglial activation and production of inflammatory mediators in neonatal rat brain following hypoxia: role of TLR4 in hypoxic microglia. *J. Neuroinflammation* 10: 23.
- Lin, X.W., et al. 2013. WW domain containing E3 ubiquitin protein ligase 1 (WWP1) negatively regulates TLR4-mediated TNF- α and IL-6 production by proteasomal degradation of TNF receptor associated factor 6 (TRAF6). *PLoS ONE* 8: e67633.