

ARL4D (N-14): sc-107423

BACKGROUND

ADP-ribosylation factors (ARFs) are highly conserved guanine nucleotide-binding proteins that enhance the ADP-ribosyltransferase activity of cholera toxin. ARFs are important in eukaryotic vesicular trafficking pathways and activating phospholipase D. ARL4D (ADP-ribosylation factor-like 4D), also known as ARL6 or ARF4L, is a 201 amino acid nuclear protein that is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL4D may play a role in membrane-associated intracellular trafficking and may promote ARF6 activation and modulate Actin remodeling by regulating ARNO. It is suggested that mutations of ARL4D is associated with Bardet-Biedl syndrome.

REFERENCES

1. Smith, S.A., et al. 1995. Isolation and mapping of a gene encoding a novel human ADP-ribosylation factor on chromosome 17q12-q21. *Genomics* 28: 113-115.
2. Katayama, T., et al. 1998. Expression of an ADP-ribosylation factor like gene, ARF4L, is induced after transient forebrain ischemia in the gerbil. *Brain Res. Mol. Brain Res.* 56: 66-75.
3. Lin, C.Y., et al. 2000. ARL4, an ARF-like protein that is developmentally regulated and localized to nuclei and nucleoli. *J. Biol. Chem.* 275: 37815-37823.
4. Nonaka, Y., et al. 2002. Recognition of ADP-ribosylation factor 4-like by HLA-A2-restricted and tumor-reactive cytotoxic T lymphocytes from patients with brain tumors. *Tissue Antigens* 60: 319-327.
5. Katayama, T., et al. 2004. Role of ARF4L in recycling between endosomes and the plasma membrane. *Cell. Mol. Neurobiol.* 24: 137-147.
6. Hofmann, I., et al. 2007. The Arl4 family of small G proteins can recruit the cytosolin ARF6 exchange factors to the plasma membrane. *Curr. Biol.* 17: 711-716.
7. Li, C.C., et al. 2007. ARL4D recruits cytosolin-2/ARNO to modulate actin remodeling. *Mol. Biol. Cell* 18: 4420-4437.
8. Alme, M.N., et al. 2007. Chronic fluoxetine treatment induces brain region-specific upregulation of genes associated with BDNF-induced long-term potentiation. *Neural Plast.* 2007: 26496.
9. Chi, J.H., et al. 2008. Increased expression of the glioma-associated antigen ARF4L after loss of the tumor suppressor PTEN. *Laboratory investigation. J. Neurosurg.* 108: 299-303.

CHROMOSOMAL LOCATION

Genetic locus: ARL4D (human) mapping to 17q21.31; Arl4d (mouse) mapping to 11 D.

SOURCE

ARL4D (N-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of ARL4D of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-107423 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ARL4D (N-14) is recommended for detection of ARL4D of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

ARL4D (N-14) is also recommended for detection of ARL4D in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for ARL4D siRNA (h): sc-94138, ARL4D siRNA (m): sc-141243, ARL4D shRNA Plasmid (h): sc-94138-SH, ARL4D shRNA Plasmid (m): sc-141243-SH, ARL4D shRNA (h) Lentiviral Particles: sc-94138-V and ARL4D shRNA (m) Lentiviral Particles: sc-141243-V.

Molecular Weight of ARL4D: 20 kDa.

Positive Controls: T-47D cell lysate: sc-2293 or Y79 cell lysate: sc-2240.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **ARL4D (F-2): sc-271273** or **ARL4D (H-2): sc-271274**, our highly recommended monoclonal alternatives to ARL4D (N-14).