



MAS1L (N-12): sc-107719

BACKGROUND

The proto-oncogene MAS1 is a G protein-coupled receptor located on the plasma membrane. In transfected NIH/3T3 cells, MAS1 has a weak focus-inducing activity. MAS1 is an antagonist of the AT₁ receptor (Angiotensin II type 1 receptor), inhibiting the actions of Angiotensin II. MAS1L (MAS1-like), also known as MAS-L or MRG (Mas-related G-protein coupled receptor MRG), is a 378 amino acid multi-pass membrane protein that belongs to the G-protein coupled receptor 1 (GPCR) family. MAS1 related proteins are thought to regulate nociceptor function and development, including the sensation or modulation of pain.

REFERENCES

1. Monnot, C., Weber, V., Stinnakre, J., Bihoreau, C., Teutsch, B., Corvol, P. and Clauser, E. 1991. Cloning and functional characterization of a novel Mas-related gene, modulating intracellular Angiotensin II actions. *Mol. Endocrinol.* 5: 1477-1487.
2. Dong, X., Han, S., Zylka, M.J., Simon, M.I. and Anderson, D.J. 2001. A diverse family of GPCRs expressed in specific subsets of nociceptive sensory neurons. *Cell* 106: 619-632.
3. Gustafson, E.L., Maguire, M., Campanella, M., Tarozzo, G., Jia, Y., Dong, X.W., Laverty, M., Murgolo, N., Priestley, T., Reggiani, A., Monsma, F. and Beltramo, M. 2005. Regulation of two rat Mas-related genes in a model of neuropathic pain. *Brain Res. Mol. Brain Res.* 142: 58-64.
4. Tallant, E.A., Ferrario, C.M. and Gallagher, P.E. 2005. Angiotensin-(1-7) inhibits growth of cardiac myocytes through activation of the Mas receptor. *Am. J. Physiol. Heart Circ. Physiol.* 289: H1560-H1566.
5. Tatemoto, K., Nozaki, Y., Tsuda, R., Konno, S., Tomura, K., Furuno, M., Ogasawara, H., Edamura, K., Takagi, H., Iwamura, H., Noguchi, M. and Naito, T. 2006. Immunoglobulin E-independent activation of mast cell is mediated by MRG receptors. *Biochem. Biophys. Res. Commun.* 349: 1322-1328.
6. Zhang, R., Yan, P.K., Zhou, C.H., Liao, J.Y. and Wang, M.W. 2007. Development of a homogeneous calcium mobilization assay for high throughput screening of Mas-related gene receptor agonists. *Acta Pharmacol. Sin.* 28: 125-131.
7. Liu, Y., Yang, F.C., Okuda, T., Dong, X., Zylka, M.J., Chen, C.L., Anderson, D.J., Kuner, R. and Ma, Q. 2008. Mechanisms of compartmentalized expression of MRG class G-protein-coupled sensory receptors. *J. Neurosci.* 28: 125-132.
8. Cox, P.J., Pitcher, T., Trim, S.A., Bell, C.H., Qin, W. and Kinloch, R.A. 2008. The effect of deletion of the orphan G-protein coupled receptor (GPCR) gene MRGE on pain-like behaviours in mice. *Mol. Pain* 4: 2.
9. Hager, U.A., Hein, A., Lennerz, J.K., Zimmermann, K., Neuhuber, W.L. and Reeh, P.W. 2008. Morphological characterization of rat Mas-related G-protein-coupled receptor C and functional analysis of agonists. *Neuroscience* 151: 242-254.

CHROMOSOMAL LOCATION

Genetic locus: MAS1L (human) mapping to 6p22.1.

SOURCE

MAS1L (N-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an N-terminal extracellular domain of MAS1L of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-107719 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MAS1L (N-12) is recommended for detection of MAS1L of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for MAS1L siRNA (h): sc-95358, MAS1L shRNA Plasmid (h): sc-95358-SH and MAS1L shRNA (h) Lentiviral Particles: sc-95358-V.

Molecular Weight of MAS1L: 42 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.