

# DAZAP2 (P-12): sc-109250

## BACKGROUND

DAZAP2 (DAZ-associated protein 2), also known as deleted in azoospermia-associated protein 2 or PRTB, is a 108 amino acid proline-rich protein that interacts with DAZ (deleted in azoospermia), a gene with multiple protein products that are deleted in infertile men. Involved in spermatogenesis, RNA splicing, transcription regulation and cell signaling, DAZAP2 is known to interact with Sox-6 and DAZL and participates in the pathogenesis of multiple myeloma. A cytoplasmic protein, DAZAP2 is encoded by a gene that maps to human chromosome 12, which encodes over 1,100 genes and comprises approximately 4.5% of the human genome. Chromosome 12 is associated with a variety of diseases and afflictions, including hypochondrogenesis, achondrogenesis, Kniest dysplasia, Noonan syndrome and trisomy 12p, which causes facial developmental defects and seizure disorders.

## REFERENCES

1. Tsui, S., et al. 2000. Identification of two novel proteins that interact with germ-cell-specific RNA-binding proteins DAZ and DAZL1. *Genomics* 65: 266-273.
2. Delgado Carrasco, J., et al. 2001. Achondrogenesis type II-hypochondrogenesis: radiological features. Case report. *An. Esp. Pediatr.* 55: 553-557.
3. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607431. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Yokoyama, T., et al. 2003. A case of Kniest dysplasia with retinal detachment and the mutation analysis. *Am. J. Ophthalmol.* 136: 1186-1188.
5. Cohen-Barak, O., et al. 2003. Sox6 regulation of cardiac myocyte development. *Nucleic Acids Res.* 31: 5941-5948.
6. Shi, Y., et al. 2004. The structure, expression and function prediction of DAZAP2, a down-regulated gene in multiple myeloma. *Genomics Proteomics Bioinformatics* 2: 47-54.
7. Forzano, F., et al. 2007. A familial case of achondrogenesis type II caused by a dominant COL2A1 mutation and "patchy" expression in the mosaic father. *Am. J. Med. Genet. A* 143A: 2815-2820.
8. Shi, Y.W., et al. 2007. Molecular features and expression of DAZAP2 in human multiple myeloma. *Chin. Med. J.* 120: 1659-1665.
9. Kim, J.E., et al. 2008. Proline-rich transcript in brain protein induces stress granule formation. *Mol. Cell. Biol.* 28: 803-813.

## CHROMOSOMAL LOCATION

Genetic locus: DAZAP2 (human) mapping to 12q13.13; Dazap2 (mouse) mapping to 15 F1.

## SOURCE

DAZAP2 (P-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of DAZAP2 of human origin.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-109250 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS


DAZAP2 (P-12) is recommended for detection of DAZAP2 of mouse and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for DAZAP2 siRNA (h): sc-95755, DAZAP2 siRNA (m): sc-142878, DAZAP2 shRNA Plasmid (h): sc-95755-SH, DAZAP2 shRNA Plasmid (m): sc-142878-SH, DAZAP2 shRNA (h) Lentiviral Particles: sc-95755-V and DAZAP2 shRNA (m) Lentiviral Particles: sc-142878-V.

Molecular Weight of DAZAP2: 17 kDa.

Positive Controls: NTERA-2 cl.D1 whole cell lysate: sc-364181.

## DATA



DAZAP2 (P-12): sc-109250. Western blot analysis of DAZAP2 expression in NTERA-2 cl.D1 whole cell lysate.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.


 MONOS  
Satisfaction  
Guaranteed

Try **DAZAP2 (G-4): sc-515182**, our highly recommended monoclonal alternative to DAZAP2 (P-12).