

MAN1A1 (D-12): sc-109313

BACKGROUND

The α -mannosidases (designated MAN1A1, MAN1A2, MAN2A1 and MAN2A2) comprise a group of soluble proteins that localize to the endoplasmic reticulum, the Golgi apparatus or the cytoplasm. Depending on their cellular location, these proteins are involved in either the processing or the degradation of newly synthesized N-glycans. MAN1A1 (mannosidase, α , class 1A, member 1), also known as MAN9, HUMM3 or HUMM9, is a 653 amino acid single-pass type II membrane protein that localizes to the Golgi apparatus and is involved in protein glycosylation. Using calcium as a cofactor, MAN1A1 functions to catalyze the hydrolysis of terminal α -D-mannose residues in select oligo-mannose oligosaccharides, a reaction that is important for the maturation of Asn-link oligosaccharides.

REFERENCES

- Blough, H.A., Pauwels, R., De Clercq, E., Cogniaux, J., Sprecher-Goldberger, S. and Thiry, L. 1986. Glycosylation inhibitors block the expression of LAV/HTLV-III (HIV) glycoproteins. *Biochem. Biophys. Res. Commun.* 141: 33-38.
- Bause, E., Bieberich, E., Rolfs, A., Völker, C. and Schmidt, B. 1993. Molecular cloning and primary structure of Man9-mannosidase from human kidney. *Eur. J. Biochem.* 217: 535-540.
- Bieberich, E. and Bause, E. 1995. Man9-mannosidase from human kidney is expressed in COS cells as a Golgi-resident type II transmembrane N-glycoprotein. *Eur. J. Biochem.* 233: 644-649.
- Tremblay, L.O., Campbell Dyke, N. and Herscovics, A. 1998. Molecular cloning, chromosomal mapping and tissue-specific expression of a novel human α 1,2-mannosidase gene involved in N-glycan maturation. *Glycobiology* 8: 585-595.
- Tremblay, L.O. and Herscovics, A. 1999. Cloning and expression of a specific human α 1,2-mannosidase that trims Man9GlcNAc2 to Man8GlcNAc2 isomer B during N-glycan biosynthesis. *Glycobiology* 9: 1073-1078.
- Gonzalez, D.S., Karaveg, K., Vandersall-Nairn, A.S., Lal, A. and Moremen, K.W. 1999. Identification, expression, and characterization of a cDNA encoding human endoplasmic reticulum mannosidase I, the enzyme that catalyzes the first mannose trimming step in mammalian Asn-linked oligosaccharide biosynthesis. *J. Biol. Chem.* 274: 21375-21386.
- Hart, M.L., Saifuddin, M. and Spear, G.T. 2003. Glycosylation inhibitors and neuraminidase enhance human immunodeficiency virus type 1 binding and neutralization by mannose-binding lectin. *J. Gen. Virol.* 84: 353-360.
- Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 604344. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: MAN1A1 (human) mapping to 6q22.31; Man1a (mouse) mapping to 10 B3.

SOURCE

MAN1A1 (D-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of MAN1A1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-109313 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MAN1A1 (D-12) is recommended for detection of MAN1A1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with family member MAN1A2.

MAN1A1 (D-12) is also recommended for detection of MAN1A1 in additional species, including porcine.

Suitable for use as control antibody for MAN1A1 siRNA (h): sc-95638, MAN1A1 siRNA (m): sc-149242, MAN1A1 shRNA Plasmid (h): sc-95638-SH, MAN1A1 shRNA Plasmid (m): sc-149242-SH, MAN1A1 shRNA (h) Lentiviral Particles: sc-95638-V and MAN1A1 shRNA (m) Lentiviral Particles: sc-149242-V.

Molecular Weight of MAN1A1: 71 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.