

POLR3B (K-14): sc-109400

BACKGROUND

RNA polymerase II (Pol II) is a multi-subunit complex responsible for catalyzing the transcription of DNA into RNA. POLR3B (polymerase (RNA) III (DNA directed) polypeptide B), also known as C128 or RPC2, is a 1,133 amino acid protein that localizes to the nucleus and belongs to the RNA polymerase β chain family. Existing as the second largest component of the multi-protein RNA polymerase III (Pol III) complex, POLR3B functions as a DNA-dependent RNA polymerase that specifically catalyzes the transcription of DNA into RNA using ribonucleoside triphosphates as substrates. The gene encoding POLR3B maps to human chromosome 12, which encodes over 1,100 genes and comprises approximately 4.5% of the human genome. Chromosome 12 is associated with a variety of diseases and afflictions, including hypochondrogenesis, achondrogenesis, Kniest dysplasia, Noonan syndrome and Trisomy 12p, which causes facial developmental defects and seizure disorders.

REFERENCES

- Jang, K.L., Collins, M.K. and Latchman, D.S. 1992. The human immunodeficiency virus Tat protein increases the transcription of human Alu repeated sequences by increasing the activity of the cellular transcription factor TFIIIC. *J. Acquir. Immune Defic. Syndr.* 5: 1142-1147.
- Allen, T.L., Brothman, A.R., Carey, J.C. and Chance, P.F. 1996. Cytogenetic and molecular analysis in Trisomy 12p. *Am. J. Med. Genet.* 63: 250-256.
- Yang, W. and Cole, W.G. 1998. Low basal transcripts of the COL2A1 collagen gene from lymphoblasts show alternative splicing of exon 12 in the Kniest form of spondyloepiphyseal dysplasia. *Hum. Mutat. Suppl.* 1: S1-S2.
- Hu, P., Wu, S., Sun, Y., Yuan, C.C., Kobayashi, R., Myers, M.P. and Hernandez, N. 2002. Characterization of human RNA polymerase III identifies orthologues for *Saccharomyces cerevisiae* RNA polymerase III subunits. *Mol. Cell. Biol.* 22: 8044-8055.
- Zumkeller, W., Volleth, M., Muschke, P., Tönnies, H., Heller, A., Liehr, T., Wieacker, P. and Stumm, M. 2004. Genotype/phenotype analysis in a patient with pure and complete Trisomy 12p. *Am. J. Med. Genet. A* 129: 261-264.
- Segel, R., Peter, I., Demmer, L.A., Cowan, J.M., Hoffman, J.D. and Bianchi, D.W. 2006. The natural history of Trisomy 12p. *Am. J. Med. Genet. A* 140: 695-703.
- Shen, E., Lei, Y., Liu, Q., Zheng, Y., Song, C., Marc, J., Wang, Y., Sun, L. and Liang, Q. 2009. Identification and characterization of INMAP, a novel interphase nucleus and mitotic apparatus protein that is involved in spindle formation and cell cycle progression. *Exp. Cell Res.* 315: 1100-1116.

CHROMOSOMAL LOCATION

Genetic locus: POLR3B (human) mapping to 12q23.3; Polr3b (mouse) mapping to 10 C1.

SOURCE

POLR3B (K-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of POLR3B of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-109400 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

POLR3B (K-14) is recommended for detection of POLR3B of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

POLR3B (K-14) is also recommended for detection of POLR3B in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for POLR3B siRNA (h): sc-95786, POLR3B siRNA (m): sc-152377, POLR3B shRNA Plasmid (h): sc-95786-SH, POLR3B shRNA Plasmid (m): sc-152377-SH, POLR3B shRNA (h) Lentiviral Particles: sc-95786-V and POLR3B shRNA (m) Lentiviral Particles: sc-152377-V.

Molecular Weight of POLR3B: 128 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **POLR3B (G-6): sc-515362**, our highly recommended monoclonal alternative to POLR3B (K-14).