

AOX1 (T-15): sc-109499

BACKGROUND

The formation of free radicals is an adverse consequence of metabolism. Free radicals endanger cells by causing oxidative damage to membranes and can lead to interruption of DNA sequences, thereby potentially resulting in carcinogenesis. As a member of the molybdo-flavoenzymes family of proteins, AOX1 (Aldehyde oxidase 1) is a 1,338 amino acid cytoplasmic protein that catalyzes the oxidation of a variety of aldehydes, leading to the production of hydrogen peroxide. Under certain conditions, AOX1 can catalyze the formation of the superoxide free radical. Defects in oxygen radical metabolism have been linked to the pathogenesis of amyotrophic lateral sclerosis (ALS), an autosomal dominant neurodegenerative disorder characterized by the death of motor neurons in the spinal cord, brain and brainstem. Significantly, AOX1 is highly expressed in the ventral horn of the spinal cord and the gene that encodes AOX1 is located in a chromosomal region that is frequently found to be implicated in ALS2. This evidence suggests that AOX1 is a candidate gene for ALS2.

REFERENCES

1. Wright, R.M., et al. 1993. cDNA cloning, characterization, and tissue-specific expression of human xanthine dehydrogenase/xanthine oxidase. *Proc. Natl. Acad. Sci. USA* 90: 10690-10694.
2. Berger, R., et al. 1995. Analysis of aldehyde oxidase and xanthine dehydrogenase/oxidase as possible candidate genes for autosomal recessive familial amyotrophic lateral sclerosis. *Somat. Cell Mol. Genet.* 21: 121-131.
3. Online Mendelian Inheritance in Man, OMIM™. 1998. Johns Hopkins University, Baltimore, MD. MIM Number: 602841. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Kurosaki, M., et al. 2004. The aldehyde oxidase gene cluster in mice and rats. Aldehyde oxidase homologue 3, a novel member of the molybdo-flavoenzyme family with selective expression in the olfactory mucosa. *J. Biol. Chem.* 279: 50482-50498.
5. Terao, M., et al. 2006. Avian and canine aldehyde oxidases. Novel insights into the biology and evolution of molybdo-flavoenzymes. *J. Biol. Chem.* 281: 19748-19761.
6. Asakawa, T., et al. 2008. Properties of 130 kDa subunit of monkey aldehyde oxidase. *Biol. Pharm. Bull.* 31: 380-385.
7. Garattini, E., et al. 2008. Mammalian aldehyde oxidases: genetics, evolution and biochemistry. *Cell. Mol. Life Sci.* 65: 1019-1048.

CHROMOSOMAL LOCATION

Genetic locus: AOX1 (human) mapping to 2q33.1; Aox1 (mouse) mapping to 1 C1.3.

SOURCE

AOX1 (T-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of AOX1 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-109499 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

AOX1 (T-15) is recommended for detection of AOX1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

AOX1 (T-15) is also recommended for detection of AOX1 in additional species, including equine, canine and porcine.

Suitable for use as control antibody for AOX1 siRNA (h): sc-94924, AOX1 siRNA (m): sc-141128, AOX1 shRNA Plasmid (h): sc-94924-SH, AOX1 shRNA Plasmid (m): sc-141128-SH, AOX1 shRNA (h) Lentiviral Particles: sc-94924-V and AOX1 shRNA (m) Lentiviral Particles: sc-141128-V.


Molecular Weight of AOX1: 150 kDa.

Positive Controls: WI-38 whole cell lysate: sc-364260, Sol8 cell lysate: sc-2249 or rat liver extract: sc-2395.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



AOX1 (T-15): sc-109499. Western blot analysis of AOX1 expression in WI-38 (A) and Sol8 (B) whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.