



## Topo III $\beta$ -2 (H-20): sc-11268

### BACKGROUND

DNA topoisomerases are nuclear enzymes that regulate the topological structure of DNA by transiently breaking and rejoining DNA strands. Although DNA topoisomerase I and DNA topoisomerase II relax both positive and negative supercoils, DNA topoisomerase III relaxes only negative supercoils. DNA topoisomerase III  $\alpha$  exists as a long and a short isoform, which are produced by alternative splicing. DNA topoisomerase III  $\alpha$ , which localizes to the nucleolus, is constitutively expressed and remains at high levels throughout the cell cycle in HL-60 cells. DNA topoisomerase III  $\beta$  exists as three isoforms, namely  $\beta$ 1,  $\beta$ 2, and  $\beta$ 3, also produced by alternative splicing. DNA topoisomerase III  $\beta$ 1 is expressed in testes, heart, and skeletal muscle, whereas  $\beta$ 2 is expressed in thymus, kidney, and pancreas.

### REFERENCES

1. D-Arpa, P., et al. 1988. cDNA cloning of human DNA topoisomerase I: catalytic activity of a 67.7 kDa carboxyl-terminal fragment. Proc. Natl. Acad. Sci. USA 85: 2543-2547.
2. Kunze, N., et al. 1991. Structure of the human type I DNA topoisomerase gene. J. Biol. Chem. 266: 9610-9616.
3. Hanai, R., et al. 1996. Human TOP3: A single-copy gene encoding DNA topoisomerase III. Proc. Natl. Acad. Sci. USA 93: 3653-3657.
4. Kawasaki, K., et al. 1997. One-megabase sequence analysis of the human immunoglobulin lambda gene locus. Genome Res. 7: 250-261.
5. Ng, S.W., et al. 1999. A new human topoisomerase III that interacts with SGS1 protein. Nucleic Acids Res. 27: 993-1000.
6. Lin, C.W., et al. 2000. Differential expression of human topoisomerase III alpha during the cell cycle progression in HL-60 leukemia cells and human peripheral blood lymphocytes. Exp. Cell. Res. 256: 225-236.

### CHROMOSOMAL LOCATION

Genetic locus: TOP3A (human) mapping to 17p11.2; Top3a (mouse) mapping to 11.

### SOURCE

Topo III $\beta$ -2 (H-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Topo III $\beta$ -2 of human origin.

### PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-11268 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### APPLICATIONS

Topo III $\beta$ -2 (H-20) is recommended for detection of DNA topoisomerase III $\beta$ -2 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.