

CRIK (G-14): sc-11627

BACKGROUND

Rho, the Ras-related small GTPase, is responsible for the regulation of actin-based cytoskeletal structures including stress fibers, focal adhesions, and the contractile ring apparatus. The Citron Rho-interacting kinase (CRIK) is a serine/threonine kinase that belongs to the myotonic dystrophy kinase family and is a known effector of Rho. CRIC can be alternatively spliced to produce two isoforms, CRIC and CRIC-short kinase (SK). CRIC is a 240 kDa protein that contains the kinase domain which is followed by the Citron sequence, and CRIC-SK is a 54 kDa protein that consists mostly of the kinase domain. Both isoforms are capable of phosphorylating exogenous substrates as well as autophosphorylation. The CRIC kinase domain is related to the Rho-associated kinase (ROK), which is a target for Rho and induces the formation of focal adhesions and stress fibers. CRIC is thought to regulate cytokinesis as it localizes to the cleavage furrow and midbody of HeLa cells during the contractile process.

REFERENCES

1. Kitagawa, M., et al. 1995. Purification and characterization of a fatty acid-activated protein kinase (PKN) from rat testis. *Biochem. J.* 310: 657-664.
2. Amano, M., et al. 1996. Identification of a putative target for Rho as the serine-threonine protein kinase N. *Science* 271: 648-650.
3. Leung, T., et al. 1996. The p160 RhoA-binding kinase ROK α is a member of a kinase family and is involved in the reorganization of the cytoskeleton. *Mol. Cell Biol.* 16: 5313-5327.
4. Di Cunto, F., et al. 1998. Citron rho-interacting kinase, a novel tissue-specific Ser/Thr kinase encompassing the Rho-Rac-binding protein Citron. *J. Biol. Chem.* 273: 29706-29711.
5. Madaule, P., et al. 1998. Role of citron kinase as a target of the small GTPase Rho in cytokinesis. *Nature* 394: 491-494.
6. Itoh, K., et al. 1999. An essential part for Rho-associated kinase in the transcellular invasion of tumor cells. *Nat. Med.* 5: 221-225.
7. Olivo, C., et al. 2000. Distinct involvement of cdc42 and RhoA GTPases in actin organization and cell shape in untransformed and Dbl oncogene transformed NIH3T3 cells. *Oncogene* 19: 1428-1436.

CHROMOSOMAL LOCATION

Genetic locus: CIT (human) mapping to 12q24.23; Cit (mouse) mapping to 5 F.

SOURCE

CRIC (G-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of CRIC of mouse origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-11627 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

CRIC (G-14) is recommended for detection of CRIC of mouse, rat and human origin by Western Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for citron siRNA (h): sc-39214, citron siRNA (m): sc-39215, citron shRNA Plasmid (h): sc-39214-SH, citron shRNA Plasmid (m): sc-39215-SH, citron shRNA (h) Lentiviral Particles: sc-39214-V and citron shRNA (m) Lentiviral Particles: sc-39215-V.

Molecular Weight (predicted) of CRIC isoforms: 231/54/177/237 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **CRIC (E-6): sc-390437** or **CRIC (C-5): sc-377449**, our highly recommended monoclonal alternatives to CRIC (G-14).