

p-PYK2 (Tyr 402): sc-11767

BACKGROUND

PYK2 (proline-rich tyrosine kinase 2), a putative member of the FAK family, exhibits 61% sequence identity with FAK within its kinase domain. Like FAK, PYK2 has been shown to be a cytoplasmic protein tyrosine kinase, which is a substrate for the intrinsic protein tyrosine kinase activity of pp60Src. PYK2 (also designated CAK β or RAFTK) is highly expressed in the central nervous system. PYK2 is rapidly phosphorylated on tyrosine residues in response to stimuli, which increases intracellular calcium levels and, in turn, activates members of the PKC family of kinases. Specifically, PYK2 is phosphorylated on Tyr 402 after stimulation with heregulin. This promotes the formation of a multiprotein complex that mediates the phosphorylation of p190 RhoGAP by Src. Activation of the PYK2 kinase leads to modulation of ion channel function and the activation of the MAPK signaling pathway. PYK2 also contains phosphorylation sites within the activation loop at Tyr 579 and Tyr 580 and on the potential GRB2-binding site at Tyr 881.

CHROMOSOMAL LOCATION

Genetic locus: PTK2B (human) mapping to 8p21.1; Ptk2b (mouse) mapping to 14 D1.

SOURCE

p-PYK2 (Tyr 402) is a goat polyclonal antibody raised against a short amino acid sequence containing Tyr 402 phosphorylated PYK2 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-11767 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

p-PYK2 (Tyr 402) is recommended for detection of Tyr 402 phosphorylated PYK2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

p-PYK2 (Tyr 402)-R is also recommended for detection of correspondingly phosphorylated PYK2 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for PYK2 siRNA (h): sc-36332, PYK2 siRNA (m): sc-36333, PYK2 shRNA Plasmid (h): sc-36332-SH, PYK2 shRNA Plasmid (m): sc-36333-SH, PYK2 shRNA (h) Lentiviral Particles: sc-36332-V and PYK2 shRNA (m) Lentiviral Particles: sc-36333-V


Molecular Weight of p-PYK2: 120 kDa.

Positive Controls: Jurkat + PMA cell lysate: sc-24718 or PYK2 (h): 293T Lysate: sc-115595.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

DATA



Western blot analysis of PYK2 phosphorylation in non-transfected: sc-117752 (A,D), untreated human PYK2 transfected: sc-115595 (B,E) and lambda protein phosphatase (sc-200312A) treated human PYK2 transfected: sc-115595 (C,F) 293T whole cell lysates. Antibodies tested include p-PYK2 (Tyr 402)-R: sc-11767-R (A,B,C) and PYK2 (N-19): sc-1514 (D,E,F).

SELECT PRODUCT CITATIONS

- Chen, Y., et al. 2008. A critical role of Rap1b in B-cell trafficking and marginal zone B-cell development. *Blood* 111: 4627-4636.
- Wang, M., et al. 2009. Functional activation of proline-rich tyrosine kinase2 (PYK2) in peripheral blood mononuclear cells from patients with systemic lupus erythematosus. *BMC Musculoskelet. Disord.* 10: 141.
- Wang, M., et al. 2009. Activation signal transduction by proline-rich tyrosine kinase 2 (PYK2) in peripheral blood mononuclear cells from patients with systemic lupus erythematosus. *Hybridoma* 28: 333-339.
- Wang, M., et al. 2009. Expression and activation of proline-rich tyrosine kinase 2 (PYK2) in peripheral blood mononuclear cells from patients with systemic lupus erythematosus. *J. Int. Med. Res.* 37: 1293-1300.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



MONOS
Satisfaction
Guaranteed

Try **p-PYK2 (13.Tyr 402): sc-293142** or **p-PYK2 (14F6): sc-81512**, our highly recommended monoclonal alternative to p-PYK2 (Tyr 402).