



## HCC-4 (N-14): sc-12264

### BACKGROUND

The human  $\beta$ -cytokine HCC-4 (also designated small inducible cytokine A16, CCL16, IL-10 inducible chemokine, chemokine CC-4, NCC-4, liver expressed chemokine (LEC), LMC, Monotactin-1, MTN-1, lymphocyte and monocyte chemoattractant, or LCC-1) displays chemotactic activity for monocytes and dendritic cells. The 11 kDa HCC-4 maps to chromosome 17q11.2, where three other CC chemokines, MIP1-1, HCC-2, and HCC-1 are clustered within a 40 kb region. HCC-4 expression, which is uniquely upregulated by IL-10, demonstrates chemotactic activity for monocytes. The ability of HCC-4 to markedly improve recognition of poorly immunogenic cells suggests that it could be an effective component of antitumor vaccines.

### REFERENCES

- Hedrick, J.A., Helms, A., Vicari, A. and Zlotnik, A. 1998. Characterization of a novel C-C chemokine, HCC-4, whose expression is increased by interleukin-10. *Blood* 91: 4242-4247.
- Shoudai, K., Hieshima, K., Fukuda, S., Iio, M., Miura, R., Imai, T., Yoshie, O. and Nomiyama, H. 1998. Isolation of cDNA encoding a novel human C-C chemokine NCC-4/LEC. *Biochim. Biophys. Acta* 1396: 273-277.
- Youn, B.S., Zhang, S., Broxmeyer, H.E., Antol, K., Fraser, M.J., Jr., Hangoc, G. and Kwon, B.S. 1998. Isolation and characterization of LMC, a novel lymphocyte and monocyte chemoattractant human C-C chemokine, with myelosuppressive activity. *Biochem. Biophys. Res. Commun.* 247: 217-222.
- Nomiyama, H., Fukuda, S., Iio, M., Tanase, S., Miura, R. and Yoshie, O. 1999. Organization of the chemokine gene cluster on human chromosome 17q11.2 containing the genes for C-C chemokine MIP1-1, HCC-2, HCC-1, LEC, and RANTES. *J. Interferon Cytokine Res.* 19: 227-234.
- Giovarelli, M., Cappello, P., Forni, G., Salcedo, T., Moore, P.A., LeFleur, D.W., Nardelli, B., Carlo, E.D., Lollini, P.L., Ruben, S., Ullrich, S., Garotta, G. and Musiani, P. 2000. Tumor rejection and immune memory elicited by locally released LEC chemokine are associated with an impressive recruitment of APCs, lymphocytes, and granulocytes. *J. Immunol.* 164: 3200-3206.

### CHROMOSOMAL LOCATION

Genetic locus: CCL16 (human) mapping to 17q11.2; Ccl16-ps (mouse) mapping to 11.

### SOURCE

HCC-4 (N-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of HCC-4 of human origin.

### PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-12264 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### APPLICATIONS

HCC-4 (N-14) is recommended for detection of HCC-4 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for HCC-4 siRNA (h): sc-39353.

Molecular Weight of HCC-4: 13.6 kDa.

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.