

PGE synthase (N-15): sc-12267

BACKGROUND

Prostaglandin E synthase (PGE synthase), also known as PIG12 and MGST1-L1, is a member of protein super family MAPEG, which consists of membrane associated proteins involved in eicosanoid and glutathione metabolism. The expression of this membrane-associated protein can be induced by the proinflammatory cytokine, IL-1 β . PGE synthase is expressed in seminal vesicles, deferent ducts, kidney, heart and spleen. The enzyme activity of PGE synthase in most organs is glutathione-dependent. PGE synthase may play a significant role in the progression of Alzheimer's disease. Human PGE synthase is localized to chromosome 9q34.11.

REFERENCES

- Ogino, N., Miyamoto, T., Yamamoto, S. and Hayaishi, O. 1977. Prostaglandin endoperoxide E isomerase from bovine vesicular gland microsomes, a glutathione-requiring enzyme. *J. Biol. Chem.* 252: 890-895.
- Tanaka, Y., Ward, S.L. and Smith, W.L. 1987. Immunochemical and kinetic evidence for two different prostaglandin H-prostaglandin E isomerases in sheep vesicular gland microsomes. *J. Biol. Chem.* 262: 1374-1381.
- Watanabe, K., Kurihara, K., Tokunaga, Y. and Hayaishi, O. 1997. Two types of microsomal prostaglandin E synthase: glutathione-dependent and -independent prostaglandin E synthases. *Biochem. Biophys. Res. Commun.* 235: 148-152.
- Jakobsson, P.J., et al. 1999. Common structural features of MAPEG—a widespread superfamily of membrane associated proteins with highly divergent functions in eicosanoid and glutathione metabolism. *Protein Sci.* 8: 689-692.
- Jakobsson, P.J., et al. 1999. Identification of human prostaglandin E synthase: a microsomal, glutathione-dependent, inducible enzyme, constituting a potential novel drug target. *Proc. Natl. Acad. Sci. USA* 96: 7220-7225.
- Forsberg, L., et al. 2000. Human glutathione dependent prostaglandin E synthase: gene structure and regulation. *FEBS Lett.* 471: 78-82.
- Satoh, K., et al. 2000. Expression of prostaglandin E synthase mRNA is induced in β -amyloid treated rat astrocytes. *Neurosci. Lett.* 283: 221-223.

CHROMOSOMAL LOCATION

Genetic locus: PTGES (human) mapping to 9q34.11; Ptges (mouse) mapping to 2 B.

SOURCE

PGE synthase (N-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of PGE synthase of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-12267 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

PGE synthase (N-15) is recommended for detection of PGE synthase of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PGE synthase (N-15) is also recommended for detection of PGE synthase in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for PGE synthase siRNA (h): sc-41642, PGE synthase siRNA (m): sc-41643, PGE synthase shRNA Plasmid (h): sc-41642-SH, PGE synthase shRNA Plasmid (m): sc-41643-SH, PGE synthase shRNA (h) Lentiviral Particles: sc-41642-V and PGE synthase shRNA (m) Lentiviral Particles: sc-41643-V.

Molecular Weight of PGE synthase: 17 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204, U-87 MG cell lysate: sc-2411 or Caki-1 cell lysate: sc-2224.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **PGE synthase (A-3): sc-166308** or **PGE synthase (B-6): sc-166309**, our highly recommended monoclonal alternatives to PGE synthase (N-15).