

Type I 5-phosphatase (C-20): sc-12305

BACKGROUND

The inositol polyphosphate 5-phosphatases selectively remove the phosphate from the 5-position of various phosphatidylinositols, which generate second messengers in response to extracellular signals. The type I inositol polyphosphate 5-phosphatase is predominantly expressed in heart, brain, skeletal muscle, and human placenta. The 5-phosphatase associates with 14-3-3 ζ , and, in turn, activates the hydrolysis of inositol 1,4,5-trisphosphate and inositol 1,3,4,5-tetrakisphosphate. The type II inositol polyphosphate 5-phosphatase encodes a full length protein and a splice variant, which is the predominant form. The 5-phosphatase is ubiquitously expressed, and hydrolyzes inositol 1,4,5-trisphosphate, inositol 1,3,4,5-tetrakisphosphate, and phosphatidylinositol 4,5-bisphosphate. Both type I and II 5-phosphatases are thought to regulate the level of intracellular calcium by acting as signal terminating enzymes.

REFERENCES

- Laxminarayan, K.M., Matzaris, M., Speed, C.J. and Mitchell, C.A. 1993. Purification and characterization of a 43-kDa membrane-associated inositol polyphosphate 5-phosphatase from human placenta. *J. Biol. Chem.* 268: 4968-4974.
- Laxminarayan, K.M., Chan, B.K., Tetaz, T., Bird, P.I. and Mitchell, C.A. 1994. Characterization of a cDNA encoding the 43-kDa membrane-associated inositol-polyphosphate 5-phosphatase. *J. Biol. Chem.* 269: 17305-17310.
- Speed, C.J., Matzaris, M., Bird, P.I. and Mitchell, C.A. 1995. Tissue distribution and intracellular localisation of the 75-kDa inositol polyphosphate 5-phosphatase. *Eur. J. Biochem.* 234: 216-224.
- Mitchell, C.A., Brown, S., Campbell, J.K., Munday, A.D. and Speed, C.J. 1996. Regulation of second messengers by the inositol polyphosphate 5-phosphatases. *Biochem. Soc. Trans.* 24: 994-1000.
- Campbell, J.K., Gurung, R., Romero, S., Speed, C.J., Andrews, R.K., Berndt, M.C. and Mitchell, C.A. 1997. Activation of the 43 kDa inositol polyphosphate 5-phosphatase by 14-3-3 ζ . *Biochemistry* 36: 15363-15370.
- Zhang, X. and Majerus, P.W. 1998. Phosphatidylinositol signalling reactions. *Semin. Cell Dev. Biol.* 9: 153-160.
- Erneux, C., Govaerts, C., Communi, D. and Pessesse, X. 1998. The diversity and possible functions of the inositol polyphosphate 5-phosphatases. *Biochim. Biophys. Acta* 1436: 185-199.
- Matzaris, M., O'Malley, C.J., Badger, A., Speed, C.J., Bird, P.I. and Mitchell, C.A. 1998. Distinct membrane and cytosolic forms of inositol polyphosphate 5-phosphatase II. Efficient membrane localization requires two discrete domains. *J. Biol. Chem.* 273: 8256-8267.

CHROMOSOMAL LOCATION

Genetic locus: INPP5A (human) mapping to 10q26.3.

SOURCE

Type I 5-phosphatase (C-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of Type I 5-phosphatase of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-12305 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Type I 5-phosphatase (C-20) is recommended for detection of Type I 5-phosphatase of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Type I 5-phosphatase (C-20) is also recommended for detection of Type I 5-phosphatase in additional species, including equine, canine and bovine.

Suitable for use as control antibody for Type I 5-phosphatase siRNA (h): sc-39095, Type I 5-phosphatase shRNA Plasmid (h): sc-39095-SH and Type I 5-phosphatase shRNA (h) Lentiviral Particles: sc-39095-V.

Molecular Weight of Type I 5-phosphatase: 43 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.