

SKIV2L2 (m): 293T Lysate: sc-123569

BACKGROUND

SKIV2L2 (superkiller viralicidic activity 2-like 2) is a 1,042 amino acid protein that belongs to the helicase family and the SKI2 subfamily. The SKIV2L2 protein localizes to nucleus, contains one helicase ATP-binding domain and one helicase C-terminal domain. SKIV2L2 is identified in the spliceosome C complex and may be involved in pre-mRNA splicing. SKIV2L2 is associated with the RNA exosome complex and involved in the 3' processing of the 7S pre-RNA to the mature 5.8S rRNA. The SKIV2L2 gene is conserved in chimpanzee, cannie, bovine, mouse, rat, chicken, zebrafish, fruit fly, mosquito, *C. elegans*, *S. pombe*, *S. cerevisiae*, *K. lactis*, *E. gossypii*, *M. grisea*, *N. crassa*, *A. thaliana* and rice, and maps to human chromosome 5q11.2. A mutation in SKIV2L2 in zebrafish causes defects in cell proliferation, suggesting that SKIV2L2 plays a role in regulating melanoblast proliferation during early stages of melanocyte regeneration.

REFERENCES

1. Andersen, J.S., Lyon, C.E., Fox, A.H., Leung, A.K., Lam, Y.W., Steen, H., Mann, M. and Lamond, A.I. 2002. Directed proteomic analysis of the human nucleolus. *Curr. Biol.* 12: 1-11.
2. Scherl, A., Couté, Y., Déon, C., Callé, A., Kindbeiter, K., Sanchez, J.C., Greco, A., Hochstrasser, D. and Diaz, J.J. 2002. Functional proteomic analysis of human nucleolus. *Mol. Biol. Cell* 13: 4100-4109.
3. Gustafson, M.P., Welcker, M., Hwang, H.C. and Clurman, B.E. 2005. Zcchc8 is a glycogen synthase kinase-3 substrate that interacts with RNA-binding proteins. *Biochem. Biophys. Res. Commun.* 338: 1359-1367.
4. Yang, C.T., Hinds, A.E., Hultman, K.A. and Johnson, S.L. 2007. Mutations in *gfpt1* and SKIV2L2 cause distinct stage-specific defects in larval melanocyte regeneration in zebrafish. *PLoS Genet.* 3: e88.
5. Iwanami, N., Okada, M., Hoa, V.Q., Seo, Y., Mitani, H., Sasaki, T., Shimizu, N., Kondoh, H., Furutani-Seiki, M. and Takahama, Y. 2009. Ethylnitrosourea-induced thymus-defective mutants identify roles of KIAA1440, TRRAP, and SKIV2L2 in teleost organ development. *Eur. J. Immunol.* 39: 2606-2616.
6. Ebmeier, C.C. and Taatjes, D.J. 2010. Activator-mediator binding regulates mediator-cofactor interactions. *Proc. Natl. Acad. Sci. USA* 107: 11283-11288.
7. Umate, P., Tuteja, N. and Tuteja, R. 2011. Genome-wide comprehensive analysis of human helicases. *Commun. Integr. Biol.* 4: 118-137.

CHROMOSOMAL LOCATION

Genetic locus: Skiv2l2 (mouse) mapping to 13 D2.2.

PRODUCT

SKIV2L2 (m): 293T Lysate represents a lysate of mouse SKIV2L2 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

SKIV2L2 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive SKIV2L2 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.