



RNF39 (Y-20): sc-130038

BACKGROUND

The RING-type zinc finger motif is present in a number of viral and eukaryotic proteins and is made of a conserved cysteine-rich domain that is able to bind two zinc atoms. Proteins that contain this conserved domain are generally involved in the ubiquitination pathway of protein degradation. RNF39 (RING finger protein 39), also known as HZFW, HZF or LIRF, is a 420 amino acid protein that localizes to the cytoplasm and contains one RING-type zinc finger and one SPRY domain. Expressed in testis, RNF39 is thought to play a role in maintaining prolonged LTP (long term-potential, or the process by which synaptic strength continues to increase following chemical stimulation). Via its ability to influence the length of the LTP response, RNF39 functions to regulate early synaptic plasticity. Multiple isoforms of RNF39 exist due to alternative splicing events.

REFERENCES

- Lepourcelet, M., Andrieux, N., Giffon, T., Pichon, L., Hampe, A., Galibert, F. and Mosser, J. 1996. Systematic sequencing of the human HLA-A/HLA-F region: establishment of a cosmid contig and identification of a new gene cluster within 37 kb of sequence. *Genomics* 37: 316-326.
- Orimo, A., Yamagishi, T., Tominaga, N., Yamauchi, Y., Hishinuma, T., Okada, K., Suzuki, M., Sato, M., Nogi, Y., Suzuki, H., Inoue, S., Yoshimura, K., Shimizu, Y. and Muramatsu, M. 2000. Molecular cloning of testis-abundant finger protein/RING finger protein 23 (RNF23), a novel RING-B box-coiled coil-B30.2 protein on the class I region of the human MHC. *Biochem. Biophys. Res. Commun.* 276: 45-51.
- Coriton, O., Lepourcelet, M., Hampe, A., Galibert, F. and Mosser, J. 2000. Transcriptional analysis of the 69-kb sequence centromeric to HLA-J: a dense and complex structure of five genes. *Mamm. Genome.* 11: 1127-1131.
- Hidaka, M., Caruana, G., Stanford, W.L., Sam, M., Correll, P.H. and Bernstein, A. 2000. Gene trapping of two novel genes, Hzf and Hhl, expressed in hematopoietic cells. *Mech. Dev.* 90: 3-15.
- Matsuo, R., Asada, A., Fujitani, K. and Inokuchi, K. 2001. LIRF, a gene induced during hippocampal long-term potentiation as an immediate-early gene, encodes a novel RING finger protein. *Biochem. Biophys. Res. Commun.* 289: 479-484.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607524. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: RNF39 (human) mapping to 6p21.33.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

SOURCE

RNF39 (Y-20) is a purified rabbit polyclonal antibody raised against RNF39 of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml PBS with < 0.1% sodium azide, 0.1% gelatin and < 0.02% sucrose.

APPLICATIONS

RNF39 (Y-20) is recommended for detection of RNF39 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for RNF39 siRNA (h): sc-95580, RNF39 shRNA Plasmid (h): sc-95580-SH and RNF39 shRNA (h) Lentiviral Particles: sc-95580-V.


Molecular Weight of RNF39: 46 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



RNF39 (Y-20): sc-130038. Western blot analysis of RNF39 expression in Hep G2 whole cell lysate.

RESEARCH USE

For research use only, not for use in diagnostic procedures.