

# SPNS1 (N-18): sc-130101

## BACKGROUND

SPNS1 (spinster homolog 1), also known as LAT, nrs or SPINL, is a 528 amino acid multi-pass membrane protein that localizes to the inner mitochondrial membrane and belongs to the spinster subfamily of the major facilitator superfamily. Existing as four alternatively spliced isoforms, SPNS1 interacts with Bcl-x and Bcl-2 and, via this interaction, is thought to be involved in necrotic or autophagic cell death. The gene encoding SPNS1 maps to human chromosome 16, which encodes over 900 genes and comprises nearly 3% of the human genome. The GAN gene is located on chromosome 16 and, with mutation, may lead to giant axonal neuropathy, a nervous system disorder characterized by increasing malfunction with growth. The rare disorder Rubinstein-Taybi syndrome is also associated with chromosome 16, as is Crohn's disease, which is a gastrointestinal inflammatory condition.

## REFERENCES

- Gilbert, F. 1999. Disease genes and chromosomes: disease maps of the human genome. *Chromosome 16. Genet. Test.* 3: 243-254.
- Nakano, Y., Fujitani, K., Kurihara, J., Ragan, J., Usui-Aoki, K., Shimoda, L., Lukacsovich, T., Suzuki, K., Sezaki, M., Sano, Y., Ueda, R., Awano, W., Kaneda, M., Umeda, M. and Yamamoto, D. 2001. Mutations in the novel membrane protein spinster interfere with programmed cell death and cause neural degeneration in *Drosophila melanogaster*. *Mol. Cell. Biol.* 21: 3775-3788.
- Yanagisawa, H., Miyashita, T., Nakano, Y. and Yamamoto, D. 2003. HSpin1, a transmembrane protein interacting with Bcl-2/Bcl-x<sub>L</sub>, induces a caspase-independent autophagic cell death. *Cell Death Differ.* 10: 798-807.
- An, H., Morrell, J.L., Jennings, J.L., Link, A.J. and Gould, K.L. 2004. Requirements of fission yeast septins for complex formation, localization, and function. *Mol. Biol. Cell* 15: 5551-5564.
- Rakha, E.A., Green, A.R., Powe, D.G., Roylance, R. and Ellis, I.O. 2006. Chromosome 16 tumor-suppressor genes in breast cancer. *Genes Chromosomes Cancer* 45: 527-535.
- Gervasini, C., Castronovo, P., Bentivegna, A., Mottadelli, F., Faravelli, F., Giovannucci-Uzielli, M.L., Pessagno, A., Lucci-Cordisco, E., Pinto, A.M., Salviati, L., Selicorni, A., Tenconi, R., Neri, G. and Larizza, L. 2007. High frequency of mosaic CREBBP deletions in Rubinstein-Taybi syndrome patients and mapping of somatic and germ-line breakpoints. *Genomics* 90: 567-573.

## CHROMOSOMAL LOCATION

Genetic locus: SPNS1 (human) mapping to 16p11.2; Spns1 (mouse) mapping to 7 F4.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.

## SOURCE

SPNS1 (N-18) is a purified rabbit polyclonal antibody raised against a peptide mapping at the N-terminus of SPNS1 of human origin.

## PRODUCT

Each vial contains 100 µg IgG in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

SPNS1 (N-18) is recommended for detection of SPNS1 of mouse and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500).

Suitable for use as control antibody for SPNS1 siRNA (h): sc-93215, SPNS1 siRNA (m): sc-153775, SPNS1 shRNA Plasmid (h): sc-93215-SH, SPNS1 shRNA Plasmid (m): sc-153775-SH, SPNS1 shRNA (h) Lentiviral Particles: sc-93215-V and SPNS1 shRNA (m) Lentiviral Particles: sc-153775-V.

Molecular Weight of SPNS1: 57 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941. 3) Immunohistochemistry: use ImmunoCruz™: sc-2051 or ABC: sc-2018 rabbit IgG Staining Systems.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.