

HIPK4 (N-21): sc-130149

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. The homeodomain-interacting protein kinases (HIPK1, HIPK2, HIPK3 and HIPK4) comprise a subfamily of kinase proteins that have a conserved protein kinase domain, as well as a separate domain that interacts with homeoproteins. HIPK4 (homeodomain-interacting protein kinase 4) is a 616 amino acid protein that localizes to the cytoplasm and exists as two isoforms that are produced by alternative splicing events. Like other members of the HIPK family, HIPK4 functions as a protein kinase that catalyzes the ATP-dependent phosphorylation of target proteins and is thought to act as a co-repressor of target transcription factors.

REFERENCES

- Kim, Y.H., Choi, C.Y., Lee, S.J., Conti, M.A. and Kim, Y. 1998. Homeodomain-interacting protein kinases, a novel family of co-repressors for homeodomain transcription factors. *J. Biol. Chem.* 273: 25875-25879.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 611712. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Isono, K., Nemoto, K., Li, Y., Takada, Y., Suzuki, R., Katsuki, M., Nakagawara, A. and Koseki, H. 2006. Overlapping roles for homeodomain-interacting protein kinases HIPK1 and HIPK2 in the mediation of cell growth in response to morphogenetic and genotoxic signals. *Mol. Cell. Biol.* 26: 2758-2771.
- Arai, S., Matsushita, A., Du, K., Yagi, K., Okazaki, Y. and Kurokawa, R. 2007. Novel homeodomain-interacting protein kinase family member, HIPK4, phosphorylates human p53 at serine 9. *FEBS Lett.* 581: 5649-5657.
- Link, N., Chen, P., Lu, W.J., Pogue, K., Chuong, A., Mata, M., Checketts, J. and Abrams, J.M. 2007. A collective form of cell death requires homeodomain interacting protein kinase. *J. Cell Biol.* 178: 567-574.
- Rinaldo, C., Siepi, F., Prodosmo, A. and Soddu, S. 2008. HIPKs: Jack of all trades in basic nuclear activities. *Biochim. Biophys. Acta.* 1783: 2124-2129.
- Boucher, M.J., Simoneau, M. and Edlund, H. 2008. The homeodomain-interacting protein kinase 2 regulates Insulin promoter factor-1/pancreatic duodenal homeobox-1 transcriptional activity. *Endocrinology* 150: 87-97.

CHROMOSOMAL LOCATION

Genetic locus: HIPK4 (human) mapping to 19q13.2.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

SOURCE

HIPK4 (N-21) is a purified rabbit polyclonal antibody raised against a peptide mapping near the N-terminus of HIPK4 of human origin.

PRODUCT

Each vial contains 100 µg of IgG in PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

HIPK4 (N-21) is recommended for detection of HIPK4 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for HIPK4 siRNA (h): sc-97189, HIPK4 shRNA Plasmid (h): sc-97189-SH and HIPK4 shRNA (h) Lentiviral Particles: sc-97189-V.


Molecular Weight of HIPK4: 80 kDa.

Positive Controls: mouse heart tissue extract.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



HIPK4 (N-21): sc-130149. Western blot analysis of HIPK4 expression in mouse heart tissue extract.

RESEARCH USE

For research use only, not for use in diagnostic procedures.