



TSSK 3 (N-22): sc-130280

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. A group of proteins that are intimately involved in this process are the serine/threonine (Ser/Thr) protein kinases. TSSK 3 (testis-specific serine/threonine-protein kinase 3), also known as SPOGA3, STK22C or STK22D, is a 268 amino acid protein that contains one protein kinase domain and belongs to the Ser/Thr protein kinase family. Using magnesium as a cofactor, TSSK 3 catalyzes the ATP-dependent phosphorylation of target proteins and is thought to be involved in male germ cell development and mature sperm function. Human TSSK 3 shares 100% homology with its mouse ortholog, implying a highly conserved function between species.

REFERENCES

- Hanks, S.K., Quinn, A.M. and Hunter, T. 1988. The protein kinase family: conserved features and deduced phylogeny of the catalytic domains. *Science* 241: 42-52.
- Hanks, S.K. and Quinn, A.M. 1991. Protein kinase catalytic domain sequence database: identification of conserved features of primary structure and classification of family members. *Meth. Enzymol.* 200: 38-62.
- Kueng, P., Nikolova, Z., Djonov, V., Hemphill, A., Rohrbach, V., Boehlen, D., Zuercher, G., Andres, A.C. and Ziemiecki, A. 1997. A novel family of serine/threonine kinases participating in spermiogenesis. *J. Cell Biol.* 139: 1851-1859.
- Zuercher, G., Rohrbach, V., Andres, A.C. and Ziemiecki, A. 2000. A novel member of the testis specific serine kinase family, TSSK 3, expressed in the Leydig cells of sexually mature mice. *Mech. Dev.* 93: 175-177.
- Visconti, P.E., Hao, Z., Purdon, M.A., Stein, P., Balsara, B.R., Testa, J.R., Herr, J.C., Moss, S.B. and Kopf, G.S. 2001. Cloning and chromosomal localization of a gene encoding a novel serine/threonine kinase belonging to the subfamily of testis-specific kinases. *Genomics* 77: 163-170.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607660. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Hao, Z., Jha, K.N., Kim, Y.H., Vemuganti, S., Westbrook, V.A., Chertihin, O., Markgraf, K., Flickinger, C.J., Coppola, M., Herr, J.C. and Visconti, P.E. 2004. Expression analysis of the human testis-specific serine/threonine kinase (TSSK) homologues. A TSSK member is present in the equatorial segment of human sperm. *Mol. Hum. Reprod.* 10: 433-444.
- Bucko-Justyna, M., Lipinski, L., Burgering, B.M. and Trzeciak, L. 2005. Characterization of testis-specific serine-threonine kinase 3 and its activation by phosphoinositide-dependent kinase-1-dependent signalling. *FEBS J.* 272: 6310-6323.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: TSSK3 (human) mapping to 1p35.1.

SOURCE

TSSK 3 (N-22) is a purified rabbit polyclonal antibody raised against a peptide mapping near the N-terminus of TSSK 3 of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

TSSK 3 (N-22) is recommended for detection of TSSK 3 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for TSSK 3 siRNA (h): sc-88779, TSSK 3 shRNA Plasmid (h): sc-88779-SH and TSSK 3 shRNA (h) Lentiviral Particles: sc-88779-V.


Molecular Weight of TSSK 3: 29 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



TSSK 3 (N-22): sc-130280. Western blot analysis of TSSK 3 expression in Jurkat whole cell lysate.

RESEARCH USE

For research use only, not for use in diagnostic procedures.