

USP21 (F-22): sc-130910

BACKGROUND

The ubiquitin (Ub) pathway involves three sequential enzymatic steps that facilitate the conjugation of Ub and Ub-like molecules to specific protein substrates. Through the use of a wide range of enzymes that can add or remove ubiquitin, the Ub pathway controls many intracellular processes such as signal transduction, transcriptional activation and cell cycle progression. USP21 (ubiquitin specific peptidase 21), also known as USP16 or USP23, is a 565 amino acid protein that belongs to the C19 peptidase family of ubiquitin carboxy-terminal hydrolases. Capable of removing ubiquitin from ubiquitinated proteins, USP21 plays a role in signal transduction and can also remove NEDD8 from NEDD8-conjugated proteins, possibly functioning to influence NEDD8-mediated protein proteolysis. Multiple isoforms of USP21 exist due to alternative splicing events.

REFERENCES

1. Yin, L., Krantz, B., Russell, N.S., Deshpande, S. and Wilkinson, K.D. 2000. Nonhydrolyzable diubiquitin analogues are inhibitors of ubiquitin conjugation and deconjugation. *Biochemistry* 39: 10001-10010.
2. Smith, T.S. and Southan, C. 2000. Sequencing, tissue distribution and chromosomal assignment of a novel ubiquitin-specific protease USP23. *Biochim. Biophys. Acta.* 1490: 184-188.
3. Gong, L., Kamitani, T., Millas, S. and Yeh, E.T. 2000. Identification of a novel isopeptidase with dual specificity for ubiquitin- and NEDD8-conjugated proteins. *J. Biol. Chem.* 275: 14212-14216.
4. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 604729. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Ismail, A. and Nawaz, Z. 2005. Nuclear hormone receptor degradation and gene transcription: an update. *IUBMB Life.* 57: 483-490.
6. Joo, H.Y., Zhai, L., Yang, C., Nie, S., Erdjument-Bromage, H., Tempst, P., Chang, C. and Wang, H. 2007. Regulation of cell cycle progression and gene expression by H2A deubiquitination. *Nature* 449: 1068-1072.
7. Nakagawa, T., Kajitani, T., Togo, S., Masuko, N., Ohdan, H., Hishikawa, Y., Koji, T., Matsuyama, T., Ikura, T., Muramatsu, M. and Ito, T. 2008. Deubiquitylation of histone H2A activates transcriptional initiation via *trans*-histone cross-talk with H3K4 di- and trimethylation. *Genes Dev.* 22: 37-49.

CHROMOSOMAL LOCATION

Genetic locus: USP21 (human) mapping to 1q23.3.

SOURCE

USP21 (F-22) is a purified rabbit polyclonal antibody raised against a peptide mapping near the C-terminus of USP21 of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

USP21 (F-22) is recommended for detection of USP21 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for USP21 siRNA (h): sc-76825, USP21 shRNA Plasmid (h): sc-76825-SH and USP21 shRNA (h) Lentiviral Particles: sc-76825-V.

Molecular Weight of USP21: 62 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941. 3) Immunohistochemistry: use ImmunoCruz™: sc-2051 or ABC: sc-2018 rabbit IgG Staining Systems.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **USP21 (3D10): sc-293400**, our highly recommended monoclonal alternative to USP21 (F-22).