

MAGE-B10 (E-14): sc-131171

BACKGROUND

The melanoma-associated antigen (MAGE) family consists of a number of antigens recognized by cytotoxic T lymphocytes. The MAGE genes were initially isolated from different kinds of tumors and, based on their virtually exclusive tumor-specific expression in adult tissues, they have been used as targets for cancer immunotherapy. MAGE genes encode for tumor-rejection antigens that are expressed in tumors of different histologic types and in normal testis and placenta. MAGE-B10 (melanoma-associated antigen B10) is a 347 amino acid protein that contains one MAGE domain. The gene encoding MAGE-B10 maps to human chromosome X, which consists of about 153 million base pairs and nearly 1,000 genes. Color blindness, hemophilia and Duchenne muscular dystrophy are well known X chromosome-linked conditions which affect males more frequently, as males carry a single X chromosome.

REFERENCES

- Chomez, P., et al. 2001. An overview of the MAGE gene family with the identification of all human members of the family. *Cancer Res.* 61: 5544-5551.
- Kirkin, A.F., et al. 2002. Cancer/testis antigens: structural and immunobiological properties. *Cancer Invest.* 20: 222-236.
- Deeb, S.S. 2005. The molecular basis of variation in human color vision. *Clin. Genet.* 67: 369-377.
- Helderman-van den Enden, A.T., et al. 2009. Recurrence risk due to germ line mosaicism: Duchenne and Becker muscular dystrophy. *Clin. Genet.* 75: 465-472.
- Kasper, C.K., et al. 2009. Mosaicism and haemophilia. *Haemophilia*. E-published.
- Online Mendelian Inheritance in Man, OMIM[™]. 2009. Johns Hopkins University, Baltimore, MD. MIM Number: 300761. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Marcar, L., et al. 2010. Mage-A cancer/testis antigens inhibit p53 function by blocking its interaction with chromatin. *Cancer Res.* 70: 10362-10370.
- Sang, M., et al. 2011. Melanoma-associated antigen genes-an update. *Cancer Lett.* 302: 85-90.
- Feng, Y., et al. 2011. When MAGE meets RING: insights into biological functions of MAGE proteins. *Protein Cell* 2: 7-12.

CHROMOSOMAL LOCATION

Genetic locus: MAGEB10 (human) mapping to Xp21.3.

SOURCE

MAGE-B10 (E-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of MAGE-B10 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-131171 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MAGE-B10 (E-14) is recommended for detection of MAGE-B10 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other MAGE family members.

Suitable for use as control antibody for MAGE-B10 siRNA (h): sc-91279, MAGE-B10 shRNA Plasmid (h): sc-91279-SH and MAGE-B10 shRNA (h) Lentiviral Particles: sc-91279-V.

Molecular Weight of MAGE-B10: 39 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.