

POLR2E (T-13): sc-131331

BACKGROUND

RNA polymerase II (Pol II) is a multi-subunit complex responsible for catalyzing the transcription of DNA into RNA. POLR2E (polymerase (RNA) II (DNA directed) polypeptide E), also designated RPB5, XAP4, RPABC1, hRPB25 or hsRPB5, is a 210 amino acid nuclear protein belonging to the archaeal rpoH/eukaryotic RPB5 RNA polymerase subunit family. POLR2E is a DNA-dependent RNA polymerase that catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. POLR2E is a component of the RNA polymerase I (Pol I), RNA polymerase II (Pol II) and RNA polymerase III (Pol III) complexes. The Pol complexes synthesize both mRNA and ribosomal RNA precursors, many functional non-coding RNAs and small RNAs including 5S rRNA and tRNAs. POLR2E interacts with TFIIF RAP 30 and is important for the association between Pol II and TFIIF.

REFERENCES

- Cheong, J.H., et al. 1995. Human RPB5, a subunit shared by eukaryotic nuclear RNA polymerases, binds human hepatitis B virus X protein and may play a role in X transactivation. *EMBO J.* 14: 143-150.
- Kershner, E., et al. 1998. Immunoaffinity purification and functional characterization of human transcription factor IIH and RNA polymerase II from clonal cell lines that conditionally express epitope-tagged subunits of the multiprotein complexes. *J. Biol. Chem.* 273: 34444-34453.
- Wei, W., et al. 2001. Direct interaction between the subunit RAP 30 of transcription factor IIF (TFIIF) and RNA polymerase subunit 5, which contributes to the association between TFIIF and RNA polymerase II. *J. Biol. Chem.* 276: 12266-12273.
- Delgermaa, L., et al. 2004. Subcellular localization of RPB5-mediating protein and its putative functional partner. *Mol. Cell. Biol.* 24: 8556-8566.
- Le, T.T., et al. 2005. Mutational analysis of human RNA polymerase II subunit 5 (RPB5): the residues critical for interactions with TFIIF subunit RAP 30 and hepatitis B virus X protein. *J. Biochem.* 138: 215-224.
- Panov, K.I., et al. 2006. RNA polymerase I-specific subunit CAST/hPAF49 has a role in the activation of transcription by upstream binding factor. *Mol. Cell. Biol.* 26: 5436-5448.
- Zaros, C., et al. 2007. Functional organization of the RPB5 subunit shared by the three yeast RNA polymerases. *Nucleic Acids Res.* 35: 634-647.
- Wada, T., et al. 2009. Screening for genetic abnormalities involved in ovarian carcinogenesis using retroviral expression libraries. *Int. J. Oncol.* 35: 973-976.
- Cloutier, P., et al. 2009. High-resolution mapping of the protein interaction network for the human transcription machinery and affinity purification of RNA polymerase II-associated complexes. *Methods* 48: 381-386.

CHROMOSOMAL LOCATION

Genetic locus: POLR2E (human) mapping to 19p13.3; Polr2e (mouse) mapping to 10 C1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

POLR2E (T-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of POLR2E of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-131331 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

POLR2E (T-13) is recommended for detection of POLR2E of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other POLR family members.

Suitable for use as control antibody for POLR2E siRNA (h): sc-97922, POLR2E siRNA (m): sc-152372, POLR2E shRNA Plasmid (h): sc-97922-SH, POLR2E shRNA Plasmid (m): sc-152372-SH, POLR2E shRNA (h) Lentiviral Particles: sc-97922-V and POLR2E shRNA (m) Lentiviral Particles: sc-152372-V.

Molecular Weight of POLR2E: 24 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **POLR2E (B-5): sc-390979** or **POLR2E (C-10): sc-390902**, our highly recommended monoclonal alternatives to POLR2E (T-13).