

# Arylsulfatase E (P-13): sc-131422

## BACKGROUND

Sulfatases hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling and degradation of macromolecules. Arylsulfatase E, also known as CDPX, CDPX1 or CDPXR, is a 589 amino acid Golgi apparatus protein expressed in the pancreas, liver and kidney. Belonging to the sulfatase family, Arylsulfatase E may be essential for the correct composition of cartilage and bone matrix during development. Inhibited by warfarin, Arylsulfatase E activity occurs under optimal pH and temperature. Defects in the gene encoding Arylsulfatase E are the cause of chondrodysplasia punctata X-linked recessive type 1 (CDPX1), which is characterized by aberrant bone mineralization, severe underdevelopment of nasal cartilage and distal phalangeal hypoplasia. Arylsulfatase E has no activity toward steroid sulfates.

## REFERENCES

- Chen, N.T., et al. 1988. Potentiation of angiostatic steroids by a synthetic inhibitor of arylsulfatase. *Lab. Invest.* 59: 453-459.
- Parenti, G., et al. 1997. The sulfatase gene family. *Curr. Opin. Genet. Dev.* 7: 386-391.
- Dahl, H.H., et al. 1999. Late diagnosis of maternal PKU in a family segregating an arylsulfatase [corrected] E mutation causing symmetrical chondrodysplasia punctata. *Mol. Genet. Metab.* 68: 503-506.
- Seidel, J., et al. 2001. Brachytelephalangi dwarfism due to the loss of ARSE and SHOX genes resulting from an X;Y translocation. *Clin. Genet.* 59: 115-121.
- Brunetti-Pierrri, N., et al. 2003. X-linked recessive chondrodysplasia punctata: spectrum of arylsulfatase E gene mutations and expanded clinical variability. *Am. J. Med. Genet. A* 117A: 164-168.
- Wolpoe, M.E., et al. 2004. Severe tracheobronchial stenosis in the X-linked recessive form of chondrodysplasia punctata. *Arch. Otolaryngol. Head Neck Surg.* 130: 1423-1426.
- Chocholska, S., et al. 2006. Molecular cytogenetic analysis of a familial interstitial deletion Xp22.2-22.3 with a highly variable phenotype in female carriers. *Am. J. Med. Genet. A* 140: 604-610.
- Nino, M., et al. 2008. Clinical and molecular analysis of arylsulfatase E in patients with brachytelephalangi chondrodysplasia punctata. *Am. J. Med. Genet. A* 146A: 997-1008.
- Casarin, A., et al. 2009. X-linked brachytelephalangi chondrodysplasia punctata: a simple trait that is not so simple. *Am. J. Med. Genet. A* 149A: 2464-2468.

## CHROMOSOMAL LOCATION

Genetic locus: ARSE (human) mapping to Xp22.33.

## SOURCE

Arylsulfatase E (P-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Arylsulfatase E of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-131422 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

Arylsulfatase E (P-13) is recommended for detection of Arylsulfatase E of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other Arylsulfatase family members .

Suitable for use as control antibody for Arylsulfatase E siRNA (h): sc-91236, Arylsulfatase E shRNA Plasmid (h): sc-91236-SH and Arylsulfatase E shRNA (h) Lentiviral Particles: sc-91236-V.

Molecular Weight of Arylsulfatase E: 66 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.