



NAP1L3 (P-14): sc-131618

BACKGROUND

Proper nucleosome assembly is critical for compacting DNA into chromatin. NAP1 (nucleosome assembly protein 1) is a nuclear protein that acts as a transcriptional regulator and functions in nucleosome assembly. NAP1L3 (nucleosome assembly protein 1-like 3), also known as MB20 or NPL3, is a 506 amino acid nuclear protein belonging to the nucleosome assembly protein (NAP) family. Expressed in human brain with weak expression in heart, NAP1L3 is encoded by a gene mapping to human chromosome Xq21.32, which is in close proximity to a region closely linked to several X-linked mental retardation syndromes. Containing nearly 153 million base pairs and housing over 1,000 genes, chromosome X acts in conjunction with chromosome Y to determine sex. There are a number of conditions related to an abnormal number and combination of sex chromosomes, some of which include Turner's syndrome, color blindness, hemophilia and Duchenne muscular dystrophy.

REFERENCES

1. Watanabe, T.K., Fujiwara, T., Nakamura, Y., Hirai, Y., Maekawa, H. and Takahashi, E. 1996. Cloning, expression pattern and mapping to Xq of NAP1L3, a gene encoding a peptide homologous to human and yeast nucleosome assembly proteins. *Cytogenet. Cell Genet.* 74: 281-285.
2. Online Mendelian Inheritance in Man, OMIM™. 1998. Johns Hopkins University, Baltimore, MD. MIM Number: 300117. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Uprichard, J. and Perry, D.J. 2002. Factor X deficiency. *Blood Rev.* 16: 97-110.
4. Sharpe, L.T., de Luca, E., Hansen, T., Jäggle, H. and Gegenfurtner, K.R. 2006. Advantages and disadvantages of human dichromacy. *J. Vis.* 6: 213-223.
5. Helderma-van den Enden, A.T., de Jong, R., den Dunnen, J.T., Houwing-Duistermaat, J.J., Kneppers, A.L., Ginjaar, H.B., Breuning, M.H. and Bakker, E. 2009. Recurrence risk due to germ line mosaicism: Duchenne and Becker muscular dystrophy. *Clin. Genet.* 75: 465-472.
6. Mullaney, R. and Murphy, D. 2009. Turner syndrome: neuroimaging findings: structural and functional. *Dev. Disabil. Res. Rev.* 15: 279-283.
7. Makishima, T., King, K., Brewer, C.C., Zalewski, C.K., Butman, J., Bakalov, V.K., Bondy, C. and Griffith, A.J. 2009. Otolaryngologic markers for the early diagnosis of Turner syndrome. *Int. J. Pediatr. Otorhinolaryngol.* 73: 1564-1567.

CHROMOSOMAL LOCATION

Genetic locus: NAP1L3 (human) mapping to Xq21.32; Nap1l3 (mouse) mapping to X E3.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SOURCE

NAP1L3 (P-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of NAP1L3 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-131618 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

NAP1L3 (P-14) is recommended for detection of NAP1L3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other NALP1L family members.

Suitable for use as control antibody for NAP1L3 siRNA (h): sc-91058, NAP1L3 siRNA (m): sc-149825, NAP1L3 shRNA Plasmid (h): sc-91058-SH, NAP1L3 shRNA Plasmid (m): sc-149825-SH, NAP1L3 shRNA (h) Lentiviral Particles: sc-91058-V and NAP1L3 shRNA (m) Lentiviral Particles: sc-149825-V.

Molecular Weight of NAP1L3: 58 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.