

USE1 (P-13): sc-132485

BACKGROUND

In eukaryotic cells, the Golgi apparatus receives newly synthesized proteins from the endoplasmic reticulum (ER) and, after covalent modification, delivers them to their destination in the cell. For membrane-directed proteins this process is believed to be carried out via vesicular transport. Correct vesicular transport is determined by specific pairing of vesicle-associated SNAREs (v-SNAREs) with those on the target membrane (t-SNAREs). Unconventional SNARE in the ER 1, also known as USE1 or protein p31, is a 259 amino acid t-SNARE that forms a larger complex with ZW10, RINT-1 and Syntaxin 18. Upon Mg²⁺-AP treatment in the presence of NSF and α -SNAP, ZW10, RINT-1 and USE1 dissociate from Syntaxin 18. USE1 is a single-pass type IV membrane protein that is localized to the endoplasmic reticulum membrane. Three named isoforms exist for USE1 as a result of alternative splicing events.

REFERENCES

- Nichols, B.J. and Pelham, H.R. 1998. SNAREs and membrane fusion in the Golgi apparatus. *Biochim. Biophys. Acta* 1404: 9-31.
- Matsuda, A., Suzuki, Y., Honda, G., Muramatsu, S., Matsuzaki, O., Nagano, Y., Doi, T., Shimotohno, K., Harada, T., Nishida, E., Hayashi, H. and Sugano, S. 2003. Large-scale identification and characterization of human genes that activate NF κ B and MAPK signaling pathways. *Oncogene* 22: 3307-3318.
- Burri, L., Varlamov, O., Doege, C.A., Hofmann, K., Beilharz, T., Rothman, J.E., Söllner, T.H. and Lithgow, T. 2003. A SNARE required for retrograde transport to the endoplasmic reticulum. *Proc. Natl. Acad. Sci. USA* 100: 9873-9877.
- Belgareh-Touzé, N., Corral-Debrinski, M., Launhardt, H., Galan, J.M., Munder, T., Le Panse, S. and Haguener-Tsapis, R. 2003. Yeast functional analysis: identification of two essential genes involved in ER to Golgi trafficking. *Traffic* 4: 607-617.
- Hirose, H., Arasaki, K., Dohmae, N., Takio, K., Hatsuzawa, K., Nagahama, M., Tani, K., Yamamoto, A., Tohyama, M. and Tagaya, M. 2004. Implication of ZW10 in membrane trafficking between the endoplasmic reticulum and Golgi. *EMBO J.* 23: 1267-1278.
- Nakajima, K., Hirose, H., Taniguchi, M., Kurashina, H., Arasaki, K., Nagahama, M., Tani, K., Yamamoto, A. and Tagaya, M. 2004. Involvement of BNIP1 in apoptosis and endoplasmic reticulum membrane fusion. *EMBO J.* 23: 3216-3226.
- Okumura, A.J., Hatsuzawa, K., Tamura, T., Nagaya, H., Saeki, K., Okumura, F., Nagao, K., Nishikawa, M., Yoshimura, A. and Wada, I. 2006. Involvement of a novel Q-SNARE, D12, in quality control of the endomembrane system. *J. Biol. Chem.* 281: 4495-4506.
- Rizo, J. and Rosenmund, C. 2008. Synaptic vesicle fusion. *Nat. Struct. Mol. Biol.* 15: 665-674.

CHROMOSOMAL LOCATION

Genetic locus: USE1 (human) mapping to 19p13.11; Use1 (mouse) mapping to 8 B3.3.

SOURCE

USE1 (P-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of USE1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-132485 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

USE1 (P-13) is recommended for detection of USE1 isoforms 1, 2 and 3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

USE1 (P-13) is also recommended for detection of USE1 isoforms 1, 2 and 3 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for USE1 siRNA (h): sc-97257, USE1 siRNA (m): sc-154939, USE1 shRNA Plasmid (h): sc-97257-SH, USE1 shRNA Plasmid (m): sc-154939-SH, USE1 shRNA (h) Lentiviral Particles: sc-97257-V and USE1 shRNA (m) Lentiviral Particles: sc-154939-V.

Molecular Weight of USE1: 29 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.