

# ANKRD42 (G-14): sc-132621

## BACKGROUND

Ankyrins are membrane adaptor molecules that play important roles in coupling integral membrane proteins to the spectrin-based cytoskeleton network. Mutations of ankyrin genes lead to severe genetic diseases such as fatal cardiac arrhythmias and hereditary spherocytosis. ANKRD42 (ankyrin repeat domain 42), also known as SARP, is a 389 amino acid protein that exists as 2 alternatively spliced isoforms that contains 9 ANK repeats. Found at highest levels in testis and sperm, ANKRD42 is ubiquitously expressed and is encoded by a gene located on human chromosome 11. Chromosome 11 houses over 1,400 genes and comprises nearly 4% of the human genome. Jervell and Lange-Nielsen syndrome, Jacobsen syndrome, Niemann-Pick disease, hereditary angioedema and Smith-Lemli-Opitz syndrome are associated with defects in genes that map to chromosome 11.

## REFERENCES

1. Davis, J.Q. and Bennett, V. 1984. Brain ankyrin. A membrane-associated protein with binding sites for spectrin, tubulin and the cytoplasmic domain of the erythrocyte anion channel. *J. Biol. Chem.* 259: 13550-13559.
2. Bennett, V., Baines, A.J. and Davis, J.Q. 1985. Ankyrin and synapsin: spectrin-binding proteins associated with brain membranes. *J. Cell. Biochem.* 29: 157-169.
3. Hryniewicz-Jankowska, A., Czogalla, A., Bok, E. and Sikorsk, A.F. 2002. Ankyrins, multifunctional proteins involved in many cellular pathways. *Folia Histochem. Cytobiol.* 40: 239-249.
4. Hopitzan, A.A., Baines, A.J. and Kordeli, E. 2006. Molecular evolution of ankyrin: gain of function in vertebrates by acquisition of an obscurin/titin-binding-related domain. *Mol. Biol. Evol.* 23: 46-55.
5. Cai, X. and Zhang, Y. 2006. Molecular evolution of the ankyrin gene family. *Mol. Biol. Evol.* 23: 550-558.
6. Browne, G.J., Fardilha, M., Oxenham, S.K., Wu, W., Helps, N.R., da Cruz E Silva, O.A., Cohen, P.T. and da Cruz E Silva, E.F. 2007. SARP, a new alternatively spliced protein phosphatase 1 and DNA interacting protein. *Biochem. J.* 402: 187-196.
7. Bhuiyan, Z.A., Momenah, T.S., Amin, A.S., Al-Khadra, A.S., Alders, M., Wilde, A.A. and Mannens, M.M. 2008. An intronic mutation leading to incomplete skipping of exon-2 in KCNQ1 rescues hearing in Jervell and Lange-Nielsen syndrome. *Prog. Biophys. Mol. Biol.* 98: 319-327.
8. Coldren, C.D., Lai, Z., Shragg, P., Rossi, E., Glidewell, S.C., Zuffardi, O., Mattina, T., Ivy, D.D., Curfs, L.M., Mattson, S.N., Riley, E.P., Treier, M. and Grossfeld, P.D. 2009. Chromosomal microarray mapping suggests a role for BSX and neurogranin in neurocognitive and behavioral defects in the 11q terminal deletion disorder (Jacobsen syndrome). *Neurogenetics* 10: 89-95.

## CHROMOSOMAL LOCATION

Genetic locus: ANKRD42 (human) mapping to 11q14.1; Ankrd42 (mouse) mapping to 7 E1.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## SOURCE

ANKRD42 (G-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of ANKRD42 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-132621 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

ANKRD42 (G-14) is recommended for detection of ANKRD42 isoforms 1 and 2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other ANKRD family members.

ANKRD42 (G-14) is also recommended for detection of ANKRD42 isoforms 1 and 2 in additional species, including equine, canine, porcine and avian.

Suitable for use as control antibody for ANKRD42 siRNA (h): sc-96652, ANKRD42 siRNA (m): sc-141101, ANKRD42 shRNA Plasmid (h): sc-96652-SH, ANKRD42 shRNA Plasmid (m): sc-141101-SH, ANKRD42 shRNA (h) Lentiviral Particles: sc-96652-V and ANKRD42 shRNA (m) Lentiviral Particles: sc-141101-V.

Molecular Weight of ANKRD42: 43 kDa.

Positive Controls: MCF7 whole cell lysate: sc-2206.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.