

TRMT1 (E-14): sc-132857

BACKGROUND

Transfer RNA (tRNA) modifications help regulate the efficiency of mRNA translation by maintaining the correct reading frames. N(2),N(2)-dimethylguanosine tRNA methyltransferase, also known as TRMT1 or tRNA(guanine-26,N(2)-N(2)) methyltransferase, is a 659 amino acid enzyme that is responsible for tRNA modifications in eukaryotes. Using S-adenosyl-L-methionine as a methyl donor, TRMT1 dimethylates a single guanine residue at position 26 of tRNA. TRMT1, which was initially identified in yeast and *C. elegans*, has a 26% and 31% sequence identity to its yeast and *C. elegans* homologs, respectively. There are two isoforms of TRMT1 produced by alternative splicing events. The TRMT1 gene maps to chromosome 19p13.2 and mutations in this gene lead to abrogated enzyme activity and a decrease in protein levels.

REFERENCES

- Edqvist, J., Straby, K.B. and Grosjean, H. 1995. Enzymatic formation of N2,N2-dimethylguanosine in eukaryotic tRNA: importance of the tRNA architecture. *Biochimie* 77: 54-61.
- Constantinesco, F., Benachenhou, N., Motorin, Y. and Grosjean, H. 1998. The tRNA(guanine-26,N2-N2) methyltransferase (TRM1) from the hyperthermophilic archaeon *Pyrococcus furiosus*: cloning, sequencing of the gene and its expression in *Escherichia coli*. *Nucleic Acids Res.* 26: 3753-3761.
- Liu, J., Liu, J. and Straby, K.B. 1998. Point and deletion mutations eliminate one or both methyl group transfers catalysed by the yeast TRM1 encoded tRNA (m22G26)dimethyltransferase. *Nucleic Acids Res.* 26: 5102-5108.
- Björk, G.R., Durand, J.M., Hagervall, T.G., Leipuviene, R., Lundgren, H.K., Nilsson, K., Chen, P., Qian, Q. and Urbonavicius, J. 1999. Transfer RNA modification: influence on translational frameshifting and metabolism. *FEBS Lett.* 452: 47-51.
- Niederberger, C., Gräub, R., Costa, A., Desgrès, J. and Schweingruber, M.E. 1999. The tRNA N2,N2-dimethylguanosine-26 methyltransferase encoded by gene TRM1 increases efficiency of suppression of an ochre codon in *Schizosaccharomyces pombe*. *FEBS Lett.* 464: 67-70.
- Constantinesco, F., Motorin, Y. and Grosjean, H. 1999. Characterisation and enzymatic properties of tRNA(guanine 26, N (2), N (2))-dimethyltransferase (Trm1p) from *Pyrococcus furiosus*. *J. Mol. Biol.* 291: 375-392.
- Liu, J. and Straby, K.B. 2000. The human tRNA(m2)(2)G(26)dimethyltransferase: functional expression and characterization of a cloned hTRM1 gene. *Nucleic Acids Res.* 28: 3445-3451.
- Urbonavicius, J., Qian, Q., Durand, J.M., Hagervall, T.G. and Björk, G.R. 2001. Improvement of reading frame maintenance is a common function for several tRNA modifications. *EMBO J.* 20: 4863-4873.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 611669. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: TRMT1 (human) mapping to 19p13.2; Trmt1 (mouse) mapping to 8 C3.

SOURCE

TRMT1 (E-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of TRMT1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-132857 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TRMT1 (E-14) is recommended for detection of TRMT1 isoforms 1 and 2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with family members TRMT5 or TRMT6.


TRMT1 (E-14) is also recommended for detection of TRMT1 isoforms 1 and 2 in additional species, including canine, bovine and porcine.

Suitable for use as control antibody for TRMT1 siRNA (h): sc-97846, TRMT1 siRNA (m): sc-154683, TRMT1 shRNA Plasmid (h): sc-97846-SH, TRMT1 shRNA Plasmid (m): sc-154683-SH, TRMT1 shRNA (h) Lentiviral Particles: sc-97846-V and TRMT1 shRNA (m) Lentiviral Particles: sc-154683-V.

Molecular Weight of TRMT1: 72 kDa.

Positive Controls: TRMT1 (m): 293T Lysate: sc-124296 or Jurkat whole cell lysate: sc-2204.

DATA



TRMT1 (E-14): sc-132857. Western blot analysis of TRMT1 expression in non-transfected: sc-117752 (A) and mouse TRMT1 transfected: sc-124296 (B) 293T whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.