



# U11/U12 snRNP 35K (K-12): sc-132861

## BACKGROUND

Small nuclear ribonucleoproteins, also known as snRNPs, combine with other proteins to form spliceosomes, a complex that catalyzes pre-mRNA splicing. There are two types of spliceosomes: U2 and U12. The U2-type spliceosome is found in all eukaryotes and excises U2-type introns, which account for the majority of pre-mRNA introns. The U12-type spliceosome removes U12-type introns, which comprise less than 1% of all human introns. The U12-type spliceosome is comprised of the U11 and U12 snRNPs as well as the U4/U6.U5 tri-snRNP. U11 and U12 bind as a U11/U12 di-snRNP complex, which recognizes the 5' splice site of the pre-mRNA during the first steps of U12-type spliceosome formation. U11/U12 snRNPs contain several proteins, including seven that are unique to the U11/U12snRNP: 65K, 59K, 48K, 35K, 31K, 25K and 20K. U11/U12 snRNP 35K is a 246 amino acid protein localized to the nucleus, and it contains a motif known to mediate RNA binding. Two named isoforms of this protein exist as a result of alternative splicing events.

## REFERENCES

- Will, C.L., Schneider, C., MacMillan, A.M., Katopodis, N.F., Neubauer, G., Wilm, M., Lührmann, R. and Query, C.C. 2001. A novel U2 and U11/U12 snRNP protein that associates with the pre-mRNA branch site. *EMBO J.* 20: 4536-4546.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 610750: World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Will, C.L., Schneider, C., Hossbach, M., Urlaub, H., Rauhut, R., Elbashir, S., Tuschl, T. and Lührmann, R. 2004. The human 18S U11/U12 snRNP contains a set of novel proteins not found in the U2-dependent spliceosome. *RNA* 10: 929-941.
- Benecke, H., Lührmann, R. and Will, C.L. 2005. The U11/U12 snRNP 65K protein acts as a molecular bridge, binding the U12 snRNA and U11-59K protein. *EMBO J.* 24: 3057-3069.
- Golas, M.M., Sander, B., Will, C.L., Lührmann, R. and Stark, H. 2005. Major conformational change in the complex SF3b upon integration into the spliceosomal U11/U12 di-snRNP as revealed by electron cryomicroscopy. *Mol. Cell.* 17: 869-883.
- Schellenberg, M.J., Edwards, R.A., Ritchie, D.B., Kent, O.A., Golas, M.M., Stark, H., Lührmann, R., Glover, J.N. and MacMillan, A.M. 2006. Crystal structure of a core spliceosomal protein interface. *Proc. Natl. Acad. Sci. USA* 103: 1266-1271.
- Wang, H., Gao, M.X., Li, L., Wang, B., Hori, N. and Sato, K. 2007. Isolation, expression, and characterization of the human ZCRB1 gene mapped to 12q12. *Genomics* 89: 59-69.
- Turunen, J.J., Will, C.L., Grote, M., Lührmann, R. and Frilander, M.J. 2008. The U11-48K protein contacts the 5' splice site of U12-type introns and the U11-59K protein. *Mol. Cell. Biol.* 28: 3548-3560.

## CHROMOSOMAL LOCATION

Genetic locus: SNRNP35 (human) mapping to 12q24.31; Snrnp35 (mouse) mapping to 5 F.

## SOURCE

U11/U12 snRNP 35K (K-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of U11/U12 snRNP 35K of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-132861 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

U11/U12 snRNP 35K (K-12) is recommended for detection of U11/U12 snRNP 35K isoforms 1 and 2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for U11/U12 snRNP 35K siRNA (h): sc-96224, U11/U12 snRNP 35K siRNA (m): sc-154832, U11/U12 snRNP 35K shRNA Plasmid (h): sc-96224-SH, U11/U12 snRNP 35K shRNA Plasmid (m): sc-154832-SH, U11/U12 snRNP 35K shRNA (h) Lentiviral Particles: sc-96224-V and U11/U12 snRNP 35K shRNA (m) Lentiviral Particles: sc-154832-V.

Molecular Weight of U11/U12 snRNP 35K: 29 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.