



Paralemmin 2 (L-14): sc-133351

BACKGROUND

Paralemmin 2, also known as PALM2, is a member of the Paralemmin family of proteins, which also includes PalmDelphin and Paralemmin, also known as Paralemmin 1. Paralemmin 2 shares 26% amino acid identity with PalmDelphin and 37% amino acid identity with Paralemmin, a widely expressed peripheral membrane protein that is involved in cell structure and shape. Paralemmin 2 is an acidic 379 amino acid protein with a C-terminal CAAX motif and it is expressed in infantile muscle, infantile heart and human skin fibroblasts. PALM2, the gene encoding Paralemmin 2, is closely adjacent to and functionally linked to the AKAP2 gene. Through differential splicing and RNA read-through, four major protein products are produced by these two genes, namely Paralemmin 2, AKAP2, AKAP-KL and PALM2-AKAP2. Paralemmin 2 is encoded by the first eight exons.

REFERENCES

- Dong, F., Feldmesser, M., Casadevall, A. and Rubin, C.S. 1998. Molecular characterization of a cDNA that encodes six isoforms of a novel murine A kinase anchor protein. *J. Biol. Chem.* 273: 6533-6541.
- Nagase, T., Ishikawa, K., Suyama, M., Kikuno, R., Hirose, M., Miyajima, N., Tanaka, A., Kotani, H., Nomura, N. and Ohara, O. 1999. Prediction of the coding sequences of unidentified human genes. XIII. The complete sequences of 100 new cDNA clones from brain which code for large proteins *in vitro*. *DNA Res.* 6: 63-70.
- Hu, B., Copeland, N.G., Gilbert, D.J., Jenkins, N.A. and Kilimann, M.W. 2001. The Paralemmin protein family: identification of Paralemmin 2, an isoform differentially spliced to AKAP2/AKAP-KL, and of PalmDelphin, a more distant cytosolic relative. *Biochem. Biophys. Res. Commun.* 285: 1369-1376.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 604582. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Castellini, M., Wolf, L.V., Chauhan, B.K., Galileo, D.S., Kilimann, M.W., Cvekl, A. and Duncan, M.K. 2005. Palm is expressed in both developing and adult mouse lens and retina. *BMC Ophthalmol.* 5: 14.
- Zhang, J.K., Zhao, L.F., Cheng, J., Guo, J., Lun, Y.Z. and Hong, Y. 2006. Screening of genes for proteins interacting with the PS1TP5 protein of hepatitis B virus: probing a human leukocyte cDNA library using the yeast two-hybrid system. *Chin. Med. J.* 119: 1884-1891.
- Scholten, A., van Veen, T.A., Vos, M.A. and Heck, A.J. 2007. Diversity of cAMP-dependent protein kinase isoforms and their anchoring proteins in mouse ventricular tissue. *J. Proteome Res.* 6: 1705-1717.

CHROMOSOMAL LOCATION

Genetic locus: PALM2 (human) mapping to 9q31.3; Palm2 (mouse) mapping to 4 B3.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SOURCE

Paralemmin 2 (L-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Paralemmin 2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-133351 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Paralemmin 2 (L-14) is recommended for detection of Paralemmin 2 isoforms 1 and 2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with Paralemmin.

Suitable for use as control antibody for Paralemmin 2 siRNA (h): sc-92466, Paralemmin 2 siRNA (m): sc-155929, Paralemmin 2 shRNA Plasmid (h): sc-92466-SH, Paralemmin 2 shRNA Plasmid (m): sc-155929-SH, Paralemmin 2 shRNA (h) Lentiviral Particles: sc-92466-V and Paralemmin 2 shRNA (m) Lentiviral Particles: sc-155929-V.

Molecular Weight of Paralemmin 2: 35 kDa.

Positive Controls: HeLa nuclear extract: sc-2120.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.