



BACKGROUND

DPYS (dihydropyrimidinase), also known as DHPase, hydantoinase or DHP, is a 519 amino acid protein that is expressed in liver and kidney tissue and belongs to the DHOase family. Functioning as a homotetramer, DPYS uses zinc as a cofactor to catalyze the second step of reductive pyrimidine degradation, namely the conversion of 5,6-dihydrouracil to 3-ureidopropionate. DPYS is subject to post-translational carbamylation, an event which enhances the ability of DPYS to bind zinc ions. Defects in the gene encoding DPYS are the cause of DHP deficiency, an autosomal recessive disorder that is characterized by epileptic or convulsive attacks, dysmorphic features and severe developmental delay and congenital microvillous atrophy.

REFERENCES

- Berger, R., et al. 1984. Dihydropyrimidine dehydrogenase deficiency leading to thymine-uraciluria. An inborn error of pyrimidine metabolism. *Clin. Chim. Acta* 141: 227-234.
- Duran, M., et al. 1991. Dihydropyrimidinuria: a new inborn error of pyrimidine metabolism. *J. Inher. Metab. Dis.* 14: 367-370.
- Hamajima, N., et al. 1996. A novel gene family defined by human dihydropyrimidinase and three related proteins with differential tissue distribution. *Gene* 180: 157-163.
- van Gennip, A.H., et al. 1997. Dihydropyrimidinase deficiency: confirmation of the enzyme defect in dihydropyrimidinuria. *J. Inher. Metab. Dis.* 20: 339-342.
- Assmann, B., et al. 1997. Dihydropyrimidinase deficiency and congenital microvillous atrophy: coincidence or genetic relation? *J. Inher. Metab. Dis.* 20: 681-688.
- Hamajima, N., et al. 1998. Dihydropyrimidinase deficiency: structural organization, chromosomal localization, and mutation analysis of the human dihydropyrimidinase gene. *Am. J. Hum. Genet.* 63: 717-726.
- Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 222748. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: DPYS (human) mapping to 8q22.3; Dpys (mouse) mapping to 15 B3.1.

SOURCE

DPYS (Q-21) is a Protein A purified rabbit polyclonal antibody raised against synthetic DPYS peptide of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml PBS with < 0.1% sodium azide, 0.1% gelatin and < 0.02% sucrose.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

DPYS (Q-21) is recommended for detection of DPYS of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for DPYS siRNA (h): sc-77763, DPYS siRNA (m): sc-143167, DPYS shRNA Plasmid (h): sc-77763-SH, DPYS shRNA Plasmid (m): sc-143167-SH, DPYS shRNA (h) Lentiviral Particles: sc-77763-V and DPYS shRNA (m) Lentiviral Particles: sc-143167-V.

Molecular Weight of DPYS monomer: 54 kDa.


Molecular Weight of DPYS homotetramer: 215 kDa.

Positive Controls: DPYS (h): 293T Lysate: sc-114343.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



DPYS (Q-21): sc-133521. Western blot analysis of DPYS expression in non-transfected: sc-117752 (A) and human DPYS transfected: sc-114343 (B) 293T whole cell lysates.



DPYS (Q-21): sc-133521. Western blot analysis of human DPYS transfected 293T whole cell lysate.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.