

eIF4E3 (V-22): sc-133542

BACKGROUND

The initiation of protein synthesis in eukaryotic cells is regulated by interactions between protein initiation factors and RNA molecules. The eukaryotic initiation complex eIF4F exists *in vitro* as a trimeric structure composed of eIF4G, eIF4E and eIF4A. Together, these proteins allow ribosome binding to mRNA by inducing the unwinding of mRNA secondary structures. eIF4E binds to the mRNA "cap" during an early step in the initiation of protein synthesis. eIF4A acts as an ATP-dependent RNA helicase and eIF4G acts as a bridge between eIF4E, eIF4A and the eIF3 complex. eIF4E3 (eukaryotic translation initiation factor 4E type 3) is a 224 amino acid protein that has residues that are specific for interactions with eIF4G and 4E-BP1. Also, eIF4E3 can recognize and bind the 7-methylguanosine-containing mRNA cap. There are two isoforms of eIF4E3 that exist as a result of alternative splicing events.

REFERENCES

- Joshi, B., Cameron, A. and Jagus, R. 2004. Characterization of mammalian eIF4E-family members. *Eur. J. Biochem.* 271: 2189-2203.
- Hernández, G., Altmann, M., Sierra, J.M., Urlaub, H., Diez del Corral, R., Schwartz, P. and Rivera-Pomar, R. 2005. Functional analysis of seven genes encoding eight translation initiation factor 4E (eIF4E) isoforms in *Drosophila*. *Mech. Dev.* 122: 529-543.
- Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 609896. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Lin, S.L., Chang, D. and Ying, S.Y. 2007. Hyaluronan stimulates transformation of androgen-independent prostate cancer. *Carcinogenesis* 28: 310-320.
- Huang, B.P., Wang, Y., Wang, X., Wang, Z. and Proud, C.G. 2008. Blocking eukaryotic initiation factor (eIF) 4F complex formation does not inhibit the mTORC1-dependent activation of protein synthesis in cardiomyocytes. *Am. J. Physiol. Heart Circ. Physiol.* 296: H505-H514.
- Evsikov, A.V. and Marín de Esvikova, C. 2008. Evolutionary origin and phylogenetic analysis of the novel oocyte-specific eukaryotic translation initiation factor 4E in Tetrapoda. *Dev. Genes Evol.* 219: 111-118.
- Mir, M.A. and Panganiban, A.T. 2008. A protein that replaces the entire cellular eIF4F complex. *EMBO J.* 27: 3129-3139.
- Yanagiya, A., Svitkin, Y.V., Shibata, S., Mikami, S., Imataka, H. and Sonenberg, N. 2008. Requirement of RNA-binding of mammalian eIF4G1 for efficient interaction of eIF4E with mRNA cap. *Mol. Cell. Biol.* 29: 1661-1669.

CHROMOSOMAL LOCATION

Genetic locus: EIF4E3 (human) mapping to 3p13; Eif4e3 (mouse) mapping to 6 D3.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

eIF4E3 (V-22) is an affinity purified rabbit polyclonal antibody raised against synthetic eIF4E3 peptide of human origin.

PRODUCT

Each vial contains 50 µg IgG in 500 µl PBS with < 0.1% sodium azide, 0.1% gelatin and < 0.02% sucrose.

APPLICATIONS

eIF4E3 (V-22) is recommended for detection of eIF4E3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for eIF4E3 siRNA (h): sc-78455, eIF4E3 siRNA (m): sc-144620, eIF4E3 shRNA Plasmid (h): sc-78455-SH, eIF4E3 shRNA Plasmid (m): sc-144620-SH, eIF4E3 shRNA (h) Lentiviral Particles: sc-78455-V and eIF4E3 shRNA (m) Lentiviral Particles: sc-144620-V.


Molecular Weight of eIF4E3: 27 kDa.

Positive Controls: human fetal liver tissue extract.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



eIF4E3 (V-22): sc-133542. Western blot analysis of eIF4E3 expression in human fetal liver tissue extract.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.