

# MOX-1 (3-RE14): sc-134389

## BACKGROUND

Closely related homeobox proteins, MOX-1 and MOX-2, belong to a family of nonclustered, diverged homeobox genes that are expressed in overlapping patterns in the paraxial mesoderm and its derivatives. MOX-1 and MOX-2 function transiently in the formation of mesodermal and mesenchymal derivatives. Specifically, MOX-1 and MOX-2 are implicated in the early steps of mesoderm formation during gastrulation and are also involved in somatic differentiation. Significantly, MOX-1 associates more strongly with Pax-1, whereas MOX-2 preferentially associates with Pax-3. Expression of MOX-1, also known as Mesenchyme homeobox 1 and MFOX1, was first detected in the newly formed mesoderm of primitive streak stage mouse embryos. MOX-1 has been shown to be critical in axial skeleton development. The human MEOX1 gene maps to chromosome 17q21.31 and encodes the MOX-1 protein.

## REFERENCES

1. Candia, A.F., Hu, J., Crosby, J., Lalley, P.A., Noden, D., Nadeau, J.H. and Wright, C.V. 1992. MOX-1 and MOX-2 define a novel homeobox gene sub-family and are differentially expressed during early mesodermal patterning in mouse embryos. *Development* 116: 1123-1136.
2. Candia, A.F. and Wright, C.V. 1996. Differential localization of MOX-1 and MOX-2 proteins indicates distinct roles during development. *Int. J. Dev. Biol.* 40: 1179-1184.
3. Stelnicki, E.J., Komuves, L.G., Holmes, D., Clavin, W., Harrison, M.R., Adzick, N.S. and Largman, C. 1997. The human homeobox genes MSX-1, MSX-2, and MOX-1 are differentially expressed in the dermis and epidermis in fetal and adult skin. *Differentiation* 62: 33-41.
4. Online Mendelian Inheritance in Man, OMIM™. 1998. Johns Hopkins University, Baltimore, MD. MIM Number: 600147. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Mankoo, B.S., Collins, N.S., Ashby, P., Grigorieva, E., Pevny, L.H., Candia, A., Wright, C.V., Rigby, P.W. and Pachnis, V. 1999. MOX-2 is a component of the genetic hierarchy controlling limb muscle development. *Nature* 400: 69-73.
6. Stamatakis, D., Kastrinaki, M., Mankoo, B.S., Pachnis, V. and Karagogeos, D. 2001. Homeodomain proteins MOX-1 and MOX-2 associate with Pax-1 and Pax-3 transcription factors. *FEBS Lett.* 499: 274-278.
7. LocusLink Report (LocusID: 4222). <http://www.ncbi.nlm.nih.gov/LocusLink/>

## CHROMOSOMAL LOCATION

Genetic locus: MEOX1 (human) mapping to 17q21.31.

## SOURCE

MOX-1 (3-RE14) is a mouse monoclonal antibody raised against recombinant MOX-1 protein of human origin.

## PRODUCT

Each vial contains 100 µg IgG<sub>2a</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

MOX-1 (3-RE14) is recommended for detection of MOX-1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for MOX-1 siRNA (h): sc-43938, MOX-1 shRNA Plasmid (h): sc-43938-SH and MOX-1 shRNA (h) Lentiviral Particles: sc-43938-V.


Molecular Weight of MOX-1: 38 kDa.

Positive Controls: human MOX-1 transfected 293T whole cell lysate.

## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

## DATA



MOX-1 (3-RE14): sc-134389. Western blot analysis of MOX-1 expression in human MOX-1 transfected (A) and non-transfected (B) 293T whole cell lysates.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.