



MAGP-2 (M-80): sc-134438

BACKGROUND

Elastic fibers endow loose connective tissue with a resilience that complements the tensile strength of collagenous fibers. They are composed of the protein elastin and a network of 10-12 nm microfibrils, which contain several glycoproteins, including fibrillin-1, fibrillin-2 and the microfibril-associated glycoproteins MAGP-1 and MAGP-2. MAGP-2 functions in maintaining extracellular matrix homeostasis through the stabilization of type I procollagen and through the binding of fibrillins to tropoelastin in the extracellular matrix of several elastic and non-elastic tissues. MAGP-2 may function outside of its role in elastic fibers and play a role in cellular differentiation through the binding of Notch 1, which leads to the release of Notch 1 extracellular domain, the subsequent activation of its signaling pathway and the release of soluble Jagged1.

REFERENCES

- Gibson, M.A. and Cleary, E.G. 1987. The immunohistochemical localisation of microfibril-associated glycoprotein (MAGP) in elastic and non-elastic tissues. *Immunol. Cell Biol.* 65: 345-356.
- Kumaratilake, J.S., Gibson, M.A., Fanning, J.C. and Cleary, E.G. 1989. The tissue distribution of microfibrils reacting with a monospecific antibody to MAGP, the major glycoprotein antigen of elastin-associated microfibrils. *Eur. J. Cell Biol.* 50: 117-127.
- Gibson, M.A., Sandberg, L.B., Grosso, L.E. and Cleary, E.G. 1991. Complementary DNA cloning establishes microfibril-associated glycoprotein (MAGP) to be a discrete component of the elastin-associated microfibrils. *J. Biol. Chem.* 266: 7596-7601.
- Chen, Y., Faraco, J., Yin, W., Germiller, J., Francke, U. and Bonadio, J. 1993. Structure, chromosomal localization and expression pattern of the murine Magp gene. *J. Biol. Chem.* 268: 27381-27389.
- Gibson, M.A., Hatzinikolas, G., Kumaratilake, J.S., Sandberg, L.B., Nicholl, J.K., Sutherland, G.R. and Cleary, E.G. 1996. Further characterization of proteins associated with elastic fiber microfibrils including the molecular cloning of MAGP-2 (MP25). *J. Biol. Chem.* 271: 1096-1103.
- Lemaire, R., Korn, J.H., Shipley, J.M. and Lafyatis, R. 2005. Increased expression of type I collagen induced by microfibril-associated glycoprotein 2: novel mechanistic insights into the molecular basis of dermal fibrosis in scleroderma. *Arthritis. Rheum.* 52: 1812-1823.
- Nehring, L.C., Miyamoto, A., Hein, P.W., Weinmaster, G. and Shipley, J.M. 2005. The extracellular matrix protein MAGP-2 interacts with Jagged1 and induces its shedding from the cell surface. *J. Biol. Chem.* 280: 20349-20355.
- Miyamoto, A., Lau, R., Hein, P.W., Shipley, J.M. and Weinmaster, G. 2006. Microfibrillar proteins MAGP-1 and MAGP-2 induce Notch1 extracellular domain dissociation and receptor activation. *J. Biol. Chem.* 281: 10089-10097.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: MFAP5 (human) mapping to 12p13.31; Mfap5 (mouse) mapping to 6 F1.

SOURCE

MAGP-2 (M-80) is a rabbit polyclonal antibody raised against amino acids 1-80 mapping at the N-terminus of MAGP-2 of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

MAGP-2 (M-80) is recommended for detection of MAGP-2 of mouse, rat and, to a lesser extent, human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for MAGP-2 siRNA (h): sc-60982, MAGP-2 siRNA (m): sc-60983, MAGP-2 shRNA Plasmid (h): sc-60982-SH, MAGP-2 shRNA Plasmid (m): sc-60983-SH, MAGP-2 shRNA (h) Lentiviral Particles: sc-60982-V and MAGP-2 shRNA (m) Lentiviral Particles: sc-60983-V.

Molecular Weight of MAGP-2: 25 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.