

ACD (M-300): sc-134659

BACKGROUND

ACD (adrenocortical dysplasia homolog), also known as PIP1, PTOP, TPP1 or TINT1 is a 544 amino acid human homolog of the mouse Acd protein and is one of six components in the telosome/shelterin complex; a complex involved in shaping and guarding telomeres. ACD is essential for the proper assembly and stabilization of the telomere-associated complex and is able to interact directly with POT1 (protection of telomeres 1) and TRF1 (telomeric repeat binding factor 1), two additional members of the multi-protein complex. Localized to the nucleus, ACD helps to control telomere length and elongation by mediating telomerase activity and telomerase access to DNA. Through its ability to control and maintain telomere growth, ACD is thought to be involved in organogenesis. Two isoforms of ACD are expressed due to alternative splicing events.

REFERENCES

1. Ye, J.Z., et al. 2004. POT1-interacting protein PIP1: a telomere length regulator that recruits POT1 to the TIN2/TRF1 complex. *Genes Dev.* 18: 1649-1654.
2. Liu, D., et al. 2004. PTOP interacts with POT1 and regulates its localization to telomeres. *Nat. Cell Biol.* 6: 673-680.
3. de Lange, T. 2005. Shelterin: the protein complex that shapes and safeguards human telomeres. *Genes Dev.* 19: 2100-2110.
4. O'Connor, M.S., et al. 2006. A critical role for TPP1 and TIN2 interaction in high-order telomeric complex assembly. *Proc. Natl. Acad. Sci. USA* 103: 11874-11879.
5. Hockemeyer, D., et al. 2007. Telomere protection by mammalian POT1 requires interaction with TPP1. *Nat. Struct. Mol. Biol.* 14: 754-761.
6. Cristofari, G., et al. 2007. Telomerase unplugged. *ACS Chem. Biol.* 2: 155-158.
7. Wang, F., et al. 2007. The POT1-TPP1 telomere complex is a telomerase processivity factor. *Nature* 445: 506-510.
8. Xin, H., et al. 2007. TPP1 is a homologue of ciliate TEBP- β and interacts with POT1 to recruit telomerase. *Nature* 445: 559-562.

CHROMOSOMAL LOCATION

Genetic locus: Acd (mouse) mapping to 8 D3.

SOURCE

ACD (M-300) is a rabbit polyclonal antibody raised against amino acids 117-416 mapping at the C-terminus of ACD of mouse origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

ACD (M-300) is recommended for detection of ACD of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for ACD siRNA (m): sc-140802, ACD shRNA Plasmid (m): sc-140802-SH and ACD shRNA (m) Lentiviral Particles: sc-140802-V.

Molecular Weight of ACD human isoforms 1/2: 58/57 kDa.

Molecular Weight of ACD mouse isoforms 1/2/3/4/5: 45/34/24/18/8 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
 Satisfaction
 Guaranteed

Try **ACD (A-8): sc-377318**, our highly recommended monoclonal alternative to ACD (M-300).