

# SNX19 (H-300): sc-134769

## BACKGROUND

Sorting nexin (SNX) proteins are members of a large family of hydrophilic PX (phospholipid-binding motif) domain-containing proteins that interact with a variety of receptor types. SNXs are widely expressed, although the tissue distribution of each SNX mRNA varies. The ability of SNXs to bind specific phospholipids, as well as their tendency to form protein-protein complexes, suggests a role for these proteins in cellular membrane trafficking and protein sorting. SNXs may also function specifically in pro-degradative sorting, internalization, endosomal recycling or simply in endosomal sorting. SNXs partially associate with cellular membranes, despite their hydrophilic nature. The SNX19 protein is specifically involved in protein metabolism and is commonly overexpressed in hematologic neoplasms.

## REFERENCES

1. Katzmann, D.J., Odorizzi, G. and Emr, S.D. 2002. Receptor downregulation and multivesicular-body sorting. *Nat. Rev. Mol. Cell Biol.* 3: 893-905.
2. Worby, C.A. and Dixon, J.E. 2002. Sorting out the cellular functions of sorting nexins. *Nat. Rev. Mol. Cell Biol.* 3: 919-931.
3. Carlton, J.G. and Cullen, P.J. 2005. Sorting nexins. *Curr. Biol.* 15: 819-820.
4. Carlton, J., Bujny, M., Rutherford, A. and Cullen, P. 2005. Sorting nexins: unifying trends and new perspectives. *Traffic* 6: 75-82.
5. Hu, Y.F., Zhang, H.L., Cai, T., Harashima, S. and Notkins, A.L. 2005. The IA-2 interactome. *Diabetologia* 48: 2576-2581.
6. Jacques, C., Baris, O., Prunier-Mirebeau, D., Savagner, F., Rodien, P., Rohmer, V., Franc, B., Guyétant, S., Malthiery, Y. and Reynier, P. 2005. Two-step differential expression analysis reveals a new set of genes involved in thyroid oncocytic tumors. *J. Clin. Endocrinol. Metab.* 90: 2314-2320.
7. Kerr, M.C., Lindsay, M.R., Luetterforst, R., Hamilton, N., Simpson, F., Parton, R.G., Gleeson, P.A. and Teasdale, R.D. 2006. Visualisation of macropinosome maturation by the recruitment of sorting nexins. *J. Cell Sci.* 119: 3967-3980.
8. Tybäckinoja, A., Saarinen-Pihkala, U., Elonen, E. and Knuutila, S. 2006. Amplified, lost and fused genes in 11q23-25 amplicon in acute myeloid leukemia, an array-CGH study. *Genes Chromosomes Cancer* 45: 257-264.

## CHROMOSOMAL LOCATION

Genetic locus: SNX19 (human) mapping to 11q24.3; Snx19 (mouse) mapping to 9 A4.

## SOURCE

SNX19 (H-300) is a rabbit polyclonal antibody raised against amino acids 326-625 mapping within an internal region of SNX19 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## APPLICATIONS

SNX19 (H-300) is recommended for detection of SNX19 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


SNX19 (H-300) is also recommended for detection of SNX19 in additional species, including equine.

Suitable for use as control antibody for SNX19 siRNA (h): sc-61589, SNX19 siRNA (m): sc-61590, SNX19 shRNA Plasmid (h): sc-61589-SH, SNX19 shRNA Plasmid (m): sc-61590-SH, SNX19 shRNA (h) Lentiviral Particles: sc-61589-V and SNX19 shRNA (m) Lentiviral Particles: sc-61590-V.

Molecular Weight of SNX19: 109 kDa.

Positive Controls: SNX19 (h): 293T Lysate: sc-114473.

## DATA



SNX19 (H-300): sc-134769. Western blot analysis of SNX19 expression in non-transfected: sc-117752 (A) and human SNX19 transfected: sc-114473 (B) 293T whole cell lysates.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **SNX19 (B-1): sc-377404**, our highly recommended monoclonal alternative to SNX19 (H-300).