

# SCO2 (H-45): sc-134778

## BACKGROUND

The SCO1 and SCO2 protein homologs belong to the SCO1/2 family of proteins. SCO1 and SCO2 both localize to the mitochondrion and are inner membrane proteins crucial for copper insertion or transport to the active site of cytochrome c oxidase (COX). COX is a crucial component in energy production because it functions as the terminal enzyme in the respiratory chain. SCO1 is predominantly expressed in highly oxidative phosphorylation tissues such as brain, heart and muscle, while SCO2 is ubiquitously expressed. Defects in the gene encoding for SCO1 may cause cytochrome c oxidase deficiency, a heterogeneous disorder. Defects in the gene encoding for SCO2 may cause cardioencephalomyopathy with cytochrome c oxidase deficiency, a fatal infantile disorder characterized by hypertrophic cardiomyopathy, lactic acidosis and gliosis.

## REFERENCES

- Jaksch, M., Paret, C., Stucka, R., Horn, N., Müller-Höcker, J., Horvath, R., Trebesch, N., Stecker, G., Freisinger, P., Thirion, C., Müller, J., Lunckwitz, R., Rödel, G., Shoubridge, E.A. and Lochmüller, H. 2001. Cytochrome c oxidase deficiency due to mutations in SCO2, encoding a mitochondrial copper-binding protein, is rescued by copper in human myoblasts. *Hum. Mol. Genet.* 10: 3025-3035.
- Balatri, E., Banci, L., Bertini, I., Cantini, F. and Ciofi-Baffoni, S. 2003. Solution structure of SCO1: a thioredoxin-like protein involved in cytochrome c oxidase assembly. *Structure* 11: 1431-1443.
- Hornig, Y.C., Cobine, P.A., Maxfield, A.B., Carr, H.S. and Winge, D.R. 2004. Specific copper transfer from the COX17 metallochaperone to both SCO1 and COX11 in the assembly of yeast cytochrome c oxidase. *J. Biol. Chem.* 279: 35334-35340.
- Leary, S.C., Kaufman, B.A., Pellicchia, G., Guercin, G.H., Mattman, A., Jaksch, M. and Shoubridge, E.A. 2004. Human SCO1 and SCO2 have independent, cooperative functions in copper delivery to cytochrome c oxidase. *Hum. Mol. Genet.* 13: 1839-1848.
- Williams, J.C., Sue, C., Banting, G.S., Yang, H., Glerum, D.M., Hendrickson, W.A. and Schon, E.A. 2005. Crystal structure of human SCO1: implications for redox signaling by a mitochondrial cytochrome c oxidase assembly protein. *J. Biol. Chem.* 280: 15202-15211.
- Hornig, Y.C., Leary, S.C., Cobine, P.A., Young, F.B., George, G.N., Shoubridge, E.A. and Winge, D.R. 2005. Human SCO1 and SCO2 function as copper-binding proteins. *J. Biol. Chem.* 280: 34113-34122.
- Horvath, R., Freisinger, P., Rubio, R., Merl, T., Bax, R., Mayr, J.A., Shawan, Muller-Hoeker, J., Pongratz, D., Moller, L.B., Horn, N. and Jaksch, M. 2005. Congenital cataract, muscular hypotonia, developmental delay and sensorineural hearing loss associated with a defect in copper metabolism. *J. Inher. Metab. Dis.* 28: 479-492.
- Abajian, C. and Rosenzweig, A.C. 2006. Crystal structure of yeast SCO1. *J. Biol. Inorg. Chem.* 11: 459-466.
- Banci, L., Bertini, I., Calderone, V., Ciofi-Baffoni, S., Mangani, S., Martinelli, M., Palumaa, P. and Wang, S. 2006. A hint for the function of human SCO1 from different structures. *Proc. Natl. Acad. Sci. USA* 103: 8595-8600.

## CHROMOSOMAL LOCATION

Genetic locus: SCO2 (human) mapping to 22q13.33; Sco2 (mouse) mapping to 15 E3.

## SOURCE

SCO2 (H-45) is a rabbit polyclonal antibody raised against amino acids 222-266 mapping at the C-terminus of SCO2 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

SCO2 (H-45) is recommended for detection of SCO2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

SCO2 (H-45) is also recommended for detection of SCO2 in additional species, including equine, bovine and porcine.

Suitable for use as control antibody for SCO2 siRNA (h): sc-61507, SCO2 siRNA (m): sc-153264, SCO2 shRNA Plasmid (h): sc-61507-SH, SCO2 shRNA Plasmid (m): sc-153264-SH, SCO2 shRNA (h) Lentiviral Particles: sc-61507-V and SCO2 shRNA (m) Lentiviral Particles: sc-153264-V.

Molecular Weight of SCO2: 26 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.